

# STRAUMANN SIRIOS™ What's New

Version 1.0.0.1851

November, 2024



## 1. FEATURES

### 1. Installation and Updates

- **Easy Installation:** the Straumann SIRIOS™ software is designed for seamless installation on laptops meeting the minimum system requirements outlined in the eFU.
- **Automatic Updates:** users will receive software updates directly within the application, ensuring the latest features and improvements are always available.

### 2. User & Account Management

- **Administrator Control:** the software includes a default Administrator account responsible for user account management.
- **Multi-User Support:** multiple dentist accounts can be set up on the system, each with a unique login and password, enabling secure, personalized access for every user.

### 3. Patient Management

- **Confidential Patient Records:** dentists can create and manage patient profiles linked exclusively to their accounts. Patient information remains private and accessible only to the associated dentist or the Administrator.
- **Comprehensive Patient Access for Admin:** the Administrator account has full access to manage all patient data on the system, ensuring smooth administration.

### 4. Case Management

- **Customizable Case Records:** dentists can create and manage cases and orders linked to their accounts, with access limited to their specific cases and those of the Administrator.
- **Data Management:** raw scan data is stored for 10 days by default. Dentists can pin cases to extend this duration until unpinned. Cases can also be duplicated with their associated scan data.

### 5. Scan Workflows

- **Tooth-Borne Restorations:** utilize the "Mark Prep" or "Cut" tools for high-resolution margin lines on preparation scans.
- **Implant-Borne Restorations:** the "Cut" tool ensures high-resolution scans for scanbody captures.
- **ClearCorrect Integration:** supports workflows for ClearCorrect aligners, revisions and retainers.
- **coDiagnostiX Integration:** simplifies the planning process with seamless integration.
- **Quick Scan:** start a basic scan without needing a treatment plan or patient profile, ideal for streamlined scanning needs.

## 6. Enhanced Scanning Features:

- **Soft Tissue Removal:** refine scans by easily removing soft tissue from captured images.
- **Shine Control:** reduce glare and reflection for clearer scan visuals.
- **Monochrome Visualization:** switch to monochrome for a simplified view.
- **Trim Scans:** edit and adjust scans directly within the software.

## 7. Scan Review Tools

- **Monochrome Visualization:** alternate view mode for enhanced detail.
- **Occlusal Contacts Display:** review occlusal contact points with accuracy.
- **Scanbody Validation:** confirm the accuracy of scanbody captures.
- **Undercut Inspection:** identify and inspect undercuts in scans.
- **Occlusal Distance Validation:** verify occlusal space requirements.
- **Margin Line Identification:** precisely outline and confirm margin lines.

## 8. Case Export Options

- **Local Export:** save and export cases to local storage.
- **Straumann AXS Integration:** effortlessly connect and export cases to Straumann AXS for further processing.
- **Cares Connect Integration:** streamline case sharing through Cares Connect.

Order ID	Status	Due date	Patient	Summary	Submitted on	Owner	Lab
10-240809-6	Ready to send	23 Aug 2024	AKS5555, Test 1-0-0-15C	-	-	Victoria Torres (Dental Doctor of America)	-
10-240809-13	Ready to send	23 Aug 2024	Ree, Kar 1-0-0-488	-	-	Victoria Torres (Dental Doctor of America)	-
10-240806-2	Ready to send	16 Aug 2024	Dracula, Comte 1-0-0-802	-	-	Victoria Torres (Dental Doctor of America)	-
10-240809-1(1)	Ready to send	23 Aug 2024	Case, Shivani AXS 1-0-0-459	-	-	Victoria Torres (Dental Doctor of America)	-
10-240820-3	Ready to send	31 Aug 2024	event_test-archive2 1-0-0-3086	-	-	Victoria Torres (Dental Doctor of America)	-

## 9. Multi language user interface

- The IOS software is fully available on the Windows language setting of the system for the following supported languages:
  - English - EN
  - Bulgarian - BG
  - French - FR
  - Danish - DA
  - German - DE
  - Greek - EL
  - Hungarian - HU
  - Italian - IT
  - Norwegian - NO
  - Portuguese - PT (Brazilian Portuguese)
  - Spanish - ES
  - Swedish – SV
- The IOS software defaults to English if the Windows language setting of the system is not supported.

## 2. BEST PRACTICE FOR OPTIMIZED SOFTWARE PERFORMANCE

### 1. System and Memory Management

- **Daily Reboot:** restart the system once a day to clear memory and optimize performance.
- **Close Background Applications:** for optimal scanner performance, ensure all other applications are closed while scanning.

### 2. Efficient Scanning

- **Limit Scanning Duration:** aim to keep scanning time within 5–10 minutes per step for the best performance.
- **Pin Cases with Reusable Data:** for cases with valuable raw scan data, consider pinning them. This allows you to duplicate and reuse the data for alternative treatment plans without unnecessary rescanning.

### 3. Recutting Tool Usage

- **Scan Immediately After Recutting:** to avoid issues with the recutting tool, scan directly after recutting without navigating to other views first.

### 4. Security and Software Updates

- **Follow Security Protocols:** adhere to your organization's information security policies, including using a strong Windows password.
- **Stay Updated:** regularly update to the latest software release to benefit from new features, bug fixes, and security enhancements.