


STRAUMANN® BLT IMPLANT SYSTEM

DRILLING PROTOCOL




Straumann® BLT Implant Ø3.3 mm

Straumann® BLT Implant					
	Bone Level Tapered	SDI Ø2.9	BLT Ø3.3	BLT Ø4.1	BLT Ø4.8
A	Diameter	Ø2.9 mm	Ø3.3 mm	Ø4.1 mm	Ø4.8 mm
B	Platform	Ø2.8 mm	Ø3.1 mm	Ø3.7 mm	Ø4.4 mm
C	Connection diameter	SC Ø2.3 mm	NC Ø2.8 mm	RC Ø3.3 mm	
D	Connection depth	3.6 mm	4.6 mm		
E	Connection depth incl. screw hole	5.6 mm	6.9 mm		
F	Apical diameter core	Ø0.6 mm	Ø1.1 mm	Ø1.8 mm	Ø2.3 mm
G	Apical diameter threads	Ø1.5 mm	Ø2.0 mm	Ø2.7 mm	Ø3.2 mm
H	Thread pitch / flank lead / depth	0.8 mm / 20° / 0.3 mm			
L	Lengths	10 - 14 mm	8 - 18 mm		
T	Tapered part / taper	3.9 mm/11°	5.3 mm/9°		5.8 mm/9°
U	Wall thickness - top	0.6 mm		0.7 mm	1.0 mm
V	Wall thickness - mid	0.5 mm	0.6 mm		0.9 mm


Workflow for BLT Ø3.3 mm NC

Implant bed preparation, Straumann® Bone Level Tapered Ø3.3 mm, Ø4.1 mm and Ø4.8 mm.



Workflow for BLT Ø4.1 mm RC

Implant bed preparation, Straumann® Bone Level Tapered Ø4.1 mm.




Note: In soft bone and very soft bone situations with a dense cortex, it is recommended to use the Profile Drill to prepare the cortical aspect of the osteotomy.

Note: In soft bone and very soft bone situations with a dense cortex, it is recommended to use the Profile Drill to prepare the cortical aspect of the osteotomy.

IMPLANT BED PREPARATION

Implant bed preparation on the example of a Bone Level Tapered Implant $\varnothing 4.1$ mm/10 mm RC in very hard bone (type 1). For specific information about the BLT $\varnothing 2.9$ mm, please refer to section 5.1.4.


After opening the gingiva, the implant bed preparation begins with the preparation of the alveolar ridge (Step 1) and the marking of the implantation site with a Round Bur (Step 1), followed by the implant bed preparation with the X Pilot VeloDrill™ and the X VeloDrill™ (Step 2 and 3), according to the endosteal implant diameter. The implant bed is widened in the cortical layer with the BLT Profile Drill (Step 5) and the threads are precut with the BLT Tap (Step 6).



STEP 1
Prepare alveolar ridge and mark implant position

Carefully reduce and smooth a narrow tapering ridge with a large Round Bur. This will provide a flat bone surface and a sufficiently wide area of bone. Mark the implantation site determined during the implant position planning with the $\varnothing 1.4$ mm Round Bur or $\varnothing 1.6$ mm Needle Drill.

Note: This step may not be applicable or may differ depending on the clinical situation (e.g. fresh extraction socket).




STEP 2
Implant axis and depth

With the $\varnothing 2.2$ mm X Pilot VeloDrill™, mark the implant axis by drilling to a depth of about 6 mm. Insert the $\varnothing 2.2$ mm Alignment Pin to check for correct implant axis orientation.

Use the $\varnothing 2.2$ mm X Pilot VeloDrill™ to prepare the implant bed to the final preparation depth. If necessary, correct any unsatisfactory implant axis orientation.


Use the $\varnothing 2.2$ mm Alignment Pin again to check the implant axis and preparation depth.

Caution: At this point take an X-ray, particularly with vertically reduced bone availability. The Alignment Pin is inserted into the drilled area, which allows a comparative visualization of the drill hole in relation to the anatomical structures.




STEP 3
Widen implant bed to $\varnothing 2.8$ mm

Widen the implant bed with the $\varnothing 2.8$ mm X VeloDrill™. If necessary, correct the implant bed position. Use the $\varnothing 2.8$ mm Depth Gauge to check the preparation depth.




STEP 4
Widen implant bed to $\varnothing 3.5$ mm

Widen the implant bed with the $\varnothing 3.5$ mm X VeloDrill™. If necessary, correct the implant bed position. Use the $\varnothing 3.5$ mm Depth Gauge to check the preparation depth.



STEP 5
Profile drilling

Shape the coronal part of the implant bed with the $\varnothing 4.1$ mm Profile Drill with the edge of the outer rim at bone level.



STEP 6
Tap drilling


Pre-cut the threads with the $\varnothing 4.1$ mm Tap Drill over the full depth of the implant bed preparation.

Caution: Profile Drills and Taps marked with two color rings must only be used for the Bone Level Tapered Implant system.

IMPLANT PLACEMENT


A Straumann® implant can be placed with the Handpiece or manually with the Ratchet. Do not exceed the recommended maximum speed of 15 rpm when you use the Handpiece.

Note: Straumann® Bone Level Tapered Implants must be rotationally oriented for both Handpiece and Ratchet insertion (see Step 4). The following step-by-step instructions show how a Bone Level Tapered Implant is placed with the Ratchet.




STEP 1
Attach Ratchet Adapter

Hold the enclosed part of the implant carrier. Attach the Ratchet Adapter to the Loxim®. You hear a click when the Adapter is attached correctly.




STEP 2
Remove implant from the carrier

Pull down the implant carrier and, simultaneously, lift the implant out of the implant carrier (keep your arms steady).




STEP 3
Place implant

Place the implant with the Ratchet Adapter into the implant bed. Use the Ratchet to move the implant into its final position turning it clockwise.



STEP 4
Correct implant orientation

While approaching the final implant position, make sure that the height markings on the blue transfer part are oriented exactly orofacially. This positions the four protrusions of the internal connection for ideal prosthetic abutment orientation. A turn to the next marking corresponds to a vertical displacement of 0.2 mm.



STEP 5
Remove Loxim®

The Loxim® can easily be re-inserted to finish an uncompleted implant placement until the implant is fully inserted. If the implant needs to be removed during implantation surgery, the Loxim® allows for counterclockwise turns. After insertion, detach the Loxim® with the Adapter.

If an insertion torque of over 35 Ncm is achieved before the implant has assumed its final position, check that the implant bed preparation is correct to avoid bone overcompression. The Loxim® is provided with a pre-determined breaking point at 80 Ncm to prevent damage to the inner configuration of the implant, thus ensuring the integrity of the interface for mounting the prosthesis.

After breakage of the Loxim®, the remaining part of the Loxim® in the implant must be removed and the implant, if not fitted correctly, has to be unscrewed with a 48h Explantation Device. The implant bed then has to be re-prepared and a new implant inserted. For further details, please consult the brochure Guidance for Implant Removal, (USLIT.426) for USA and 152.806 (English) and 153.866 (French) for Canada.