

**System Overview**  
**BLX SURGERY**

**Implants**

Ø 3.5 mm			Ø 3.75 mm			Ø 4.0 mm			Ø 4.5 mm			Ø 5.0 mm			Ø 5.5 mm			Ø 6.5 mm					
Roxolid® SLActive®			Roxolid® SLA®			Roxolid® SLActive®			Roxolid® SLA®			Roxolid® SLActive®			Roxolid® SLA®			Roxolid® SLActive®			Roxolid® SLA®		
6 mm	–	–	6 mm	061.4306	061.4506	6 mm	061.5306	061.5506	6 mm	061.6306	061.6506	6 mm	061.7306	061.7506	6 mm	061.8306	061.8506	6 mm	061.9306	061.9506			
8 mm	061.3308	061.3508	8 mm	061.4308	061.4508	8 mm	061.5308	061.5508	8 mm	061.6308	061.6508	8 mm	061.7308	061.7508	8 mm	061.8308	061.8508	8 mm	061.9308	061.9508			
10 mm	061.3310	061.3510	10 mm	061.4310	061.4510	10 mm	061.5310	061.5510	10 mm	061.6310	061.6510	10 mm	061.7310	061.7510	10 mm	061.8310	061.8510	10 mm	061.9310	061.9510			
12 mm	061.3312	061.3512	12 mm	061.4312	061.4512	12 mm	061.5312	061.5512	12 mm	061.6312	061.6512	12 mm	061.7312	061.7512	12 mm	061.8312	061.8512	12 mm	061.9312	061.9512			
14 mm	061.3314	061.3514	14 mm	061.4314	061.4514	14 mm	061.5314	061.5514	14 mm	061.6314	061.6514	14 mm	061.7314	061.7514	14 mm	061.8314	061.8514	14 mm	061.9314	061.9514			
16 mm	061.3316	061.3516	16 mm	061.4316	061.4516	16 mm	061.5316	061.5516	16 mm	061.6316	061.6516	16 mm	061.7316	061.7516	16 mm	061.8316	061.8516	16 mm	061.9316	061.9516			
18 mm	061.3318	061.3518	18 mm	061.4318	061.4518	18 mm	061.5318	061.5518	18 mm	061.6318	061.6518	18 mm	061.7318	061.7518	18 mm	–	–	18 mm	–	–			

**Implant Drivers**



**VeloDrill™**  
**Alignment Pins**  
**Depth Gauges**

Drill No.	Pin / Gauge	Length	Part No.
Drill No. 1 (Pilot Drill)		Short, L 34 mm	066.1301
		Medium, L 38 mm	066.1501
		Long, L 42 mm	066.1701
		Alignment Pin	046.799
Drill No. 2		Short, L 34 mm	066.1302
		Medium, L 38 mm	066.1502
		Long, L 42 mm	066.1702
		Depth Gauge	046.800
Drill No. 3		Short, L 34 mm	066.1303
		Medium, L 38 mm	066.1503
		Long, L 42 mm	066.1703
		Depth Gauge	046.801
Drill No. 4		Short, L 34 mm	066.1304
		Medium, L 38 mm	066.1504
		Long, L 42 mm	066.1704
		Depth Gauge	046.802
Drill No. 5		Short, L 34 mm	066.1305
		Medium, L 38 mm	066.1505
		Long, L 42 mm	066.1705
		Depth Gauge	046.803
Drill No. 6		Short, L 34 mm	066.1306
		Medium, L 38 mm	066.1506
		Long, L 42 mm	066.1706
		Depth Gauge	046.804
Drill No. 7		Short, L 34 mm	066.1307
		Medium, L 38 mm	066.1507
		Long, L 42 mm	066.1707
		Depth Gauge	046.805
Drill No. 8		Short, L 34 mm	066.1308
		Medium, L 38 mm	066.1508
		Long, L 42 mm	–
		Depth Gauge	046.806
Drill No. 9		Short, L 34 mm	066.1309
		Medium, L 38 mm	066.1509
		Long, L 42 mm	–
		Depth Gauge	046.807

**Closure Cap**  
**H 0.4 mm**



**Healing Abutment**

Healing Abutment for Crown, AH 2 mm	Healing Abutment for Crown, AH 4 mm	Healing Abutment for Bridge, AH 2 mm	Healing Abutment for Bridge, AH 4 mm	Healing Abutments for Crown, AH 2 mm	Healing Abutments for Crown, AH 4 mm	Healing Abutments for Crown – Wide Profile, AH 2 mm	Healing Abutments for Crown – Wide Profile, AH 4 mm	Healing Abutment for Bridge, AH 2 mm	Healing Abutment for Bridge, AH 4 mm																																																																																																																																										
<table border="1"> <tr><th colspan="5">Diameter (Ø)</th></tr> <tr><th>GH</th><th>Ø 4.0 mm</th><th>Ø 5.0 mm</th><th>Ø 6.5 mm</th><th>Ø 7.5 mm</th></tr> <tr><td>1.5 mm</td><td>064.4202S</td><td>064.4212S</td><td>064.4222S</td><td>–</td></tr> <tr><td>2.5 mm</td><td>064.4204S</td><td>064.4214S</td><td>064.4224S</td><td>064.4234S*</td></tr> <tr><td>3.5 mm</td><td>064.4206S</td><td>064.4216S</td><td>064.4226S</td><td>064.4236S*</td></tr> </table>	Diameter (Ø)					GH	Ø 4.0 mm	Ø 5.0 mm	Ø 6.5 mm	Ø 7.5 mm	1.5 mm	064.4202S	064.4212S	064.4222S	–	2.5 mm	064.4204S	064.4214S	064.4224S	064.4234S*	3.5 mm	064.4206S	064.4216S	064.4226S	064.4236S*	<table border="1"> <tr><th colspan="5">Diameter (Ø)</th></tr> <tr><th>GH</th><th>Ø 4.0 mm</th><th>Ø 5.0 mm</th><th>Ø 6.5 mm</th><th>Ø 7.5 mm</th></tr> <tr><td>1.5 mm</td><td>064.4203S</td><td>064.4213S</td><td>064.4223S</td><td>–</td></tr> <tr><td>2.5 mm</td><td>064.4205S</td><td>064.4215S</td><td>064.4225S</td><td>064.4235S*</td></tr> <tr><td>3.5 mm</td><td>064.4207S</td><td>064.4217S</td><td>064.4227S</td><td>064.4237S*</td></tr> </table>	Diameter (Ø)					GH	Ø 4.0 mm	Ø 5.0 mm	Ø 6.5 mm	Ø 7.5 mm	1.5 mm	064.4203S	064.4213S	064.4223S	–	2.5 mm	064.4205S	064.4215S	064.4225S	064.4235S*	3.5 mm	064.4207S	064.4217S	064.4227S	064.4237S*	<table border="1"> <tr><th colspan="2">Diameter (Ø)</th></tr> <tr><th>GH</th><th>Ø 5.0 mm</th></tr> <tr><td>1.5 mm</td><td>064.4232S</td></tr> </table>	Diameter (Ø)		GH	Ø 5.0 mm	1.5 mm	064.4232S	<table border="1"> <tr><th colspan="2">Diameter (Ø)</th></tr> <tr><th>GH</th><th>Ø 5.0 mm</th></tr> <tr><td>1.5 mm</td><td>064.4233S</td></tr> </table>	Diameter (Ø)		GH	Ø 5.0 mm	1.5 mm	064.4233S	<table border="1"> <tr><th colspan="4">Diameter (Ø)</th></tr> <tr><th>GH</th><th>Ø 4.0 mm</th><th>Ø 5.0 mm</th><th>Ø 7.5 mm</th></tr> <tr><td>1.5 mm</td><td>064.4202S</td><td>064.4212S</td><td>064.4222S</td></tr> <tr><td>2.5 mm</td><td>064.4204S</td><td>064.4214S</td><td>064.4224S</td></tr> <tr><td>3.5 mm</td><td>064.4206S</td><td>064.4216S</td><td>064.4226S</td></tr> </table>	Diameter (Ø)				GH	Ø 4.0 mm	Ø 5.0 mm	Ø 7.5 mm	1.5 mm	064.4202S	064.4212S	064.4222S	2.5 mm	064.4204S	064.4214S	064.4224S	3.5 mm	064.4206S	064.4216S	064.4226S	<table border="1"> <tr><th colspan="4">Diameter (Ø)</th></tr> <tr><th>GH</th><th>Ø 4.0 mm</th><th>Ø 5.0 mm</th><th>Ø 7.5 mm</th></tr> <tr><td>1.5 mm</td><td>064.4203S</td><td>064.4213S</td><td>–</td></tr> <tr><td>2.5 mm</td><td>064.4205S</td><td>064.4215S</td><td>064.4225S</td></tr> <tr><td>3.5 mm</td><td>064.4207S</td><td>064.4217S</td><td>064.4227S</td></tr> </table>	Diameter (Ø)				GH	Ø 4.0 mm	Ø 5.0 mm	Ø 7.5 mm	1.5 mm	064.4203S	064.4213S	–	2.5 mm	064.4205S	064.4215S	064.4225S	3.5 mm	064.4207S	064.4217S	064.4227S	<table border="1"> <tr><th colspan="3">Diameter (Ø)</th></tr> <tr><th>GH</th><th>Ø 6.0 mm</th><th>Ø 7.5 mm</th></tr> <tr><td>0.75 mm</td><td>064.8201S</td><td>064.8203S</td></tr> <tr><td>1.5 mm</td><td>064.8212S</td><td>064.8214S</td></tr> </table>	Diameter (Ø)			GH	Ø 6.0 mm	Ø 7.5 mm	0.75 mm	064.8201S	064.8203S	1.5 mm	064.8212S	064.8214S	<table border="1"> <tr><th colspan="3">Diameter (Ø)</th></tr> <tr><th>GH</th><th>Ø 6.0 mm</th><th>Ø 7.5 mm</th></tr> <tr><td>0.75 mm</td><td>064.8202S</td><td>064.8204S</td></tr> <tr><td>1.5 mm</td><td>064.8213S</td><td>064.8215S</td></tr> </table>	Diameter (Ø)			GH	Ø 6.0 mm	Ø 7.5 mm	0.75 mm	064.8202S	064.8204S	1.5 mm	064.8213S	064.8215S	<table border="1"> <tr><th colspan="2">Diameter (Ø)</th></tr> <tr><th>GH</th><th>Ø 5.0 mm</th></tr> <tr><td>1.5 mm</td><td>064.4232S</td></tr> </table>	Diameter (Ø)		GH	Ø 5.0 mm	1.5 mm	064.4232S	<table border="1"> <tr><th colspan="2">Diameter (Ø)</th></tr> <tr><th>GH</th><th>Ø 5.0 mm</th></tr> <tr><td>1.5 mm</td><td>064.4233S</td></tr> </table>	Diameter (Ø)		GH	Ø 5.0 mm	1.5 mm	064.4233S
Diameter (Ø)																																																																																																																																																			
GH	Ø 4.0 mm	Ø 5.0 mm	Ø 6.5 mm	Ø 7.5 mm																																																																																																																																															
1.5 mm	064.4202S	064.4212S	064.4222S	–																																																																																																																																															
2.5 mm	064.4204S	064.4214S	064.4224S	064.4234S*																																																																																																																																															
3.5 mm	064.4206S	064.4216S	064.4226S	064.4236S*																																																																																																																																															
Diameter (Ø)																																																																																																																																																			
GH	Ø 4.0 mm	Ø 5.0 mm	Ø 6.5 mm	Ø 7.5 mm																																																																																																																																															
1.5 mm	064.4203S	064.4213S	064.4223S	–																																																																																																																																															
2.5 mm	064.4205S	064.4215S	064.4225S	064.4235S*																																																																																																																																															
3.5 mm	064.4207S	064.4217S	064.4227S	064.4237S*																																																																																																																																															
Diameter (Ø)																																																																																																																																																			
GH	Ø 5.0 mm																																																																																																																																																		
1.5 mm	064.4232S																																																																																																																																																		
Diameter (Ø)																																																																																																																																																			
GH	Ø 5.0 mm																																																																																																																																																		
1.5 mm	064.4233S																																																																																																																																																		
Diameter (Ø)																																																																																																																																																			
GH	Ø 4.0 mm	Ø 5.0 mm	Ø 7.5 mm																																																																																																																																																
1.5 mm	064.4202S	064.4212S	064.4222S																																																																																																																																																
2.5 mm	064.4204S	064.4214S	064.4224S																																																																																																																																																
3.5 mm	064.4206S	064.4216S	064.4226S																																																																																																																																																
Diameter (Ø)																																																																																																																																																			
GH	Ø 4.0 mm	Ø 5.0 mm	Ø 7.5 mm																																																																																																																																																
1.5 mm	064.4203S	064.4213S	–																																																																																																																																																
2.5 mm	064.4205S	064.4215S	064.4225S																																																																																																																																																
3.5 mm	064.4207S	064.4217S	064.4227S																																																																																																																																																
Diameter (Ø)																																																																																																																																																			
GH	Ø 6.0 mm	Ø 7.5 mm																																																																																																																																																	
0.75 mm	064.8201S	064.8203S																																																																																																																																																	
1.5 mm	064.8212S	064.8214S																																																																																																																																																	
Diameter (Ø)																																																																																																																																																			
GH	Ø 6.0 mm	Ø 7.5 mm																																																																																																																																																	
0.75 mm	064.8202S	064.8204S																																																																																																																																																	
1.5 mm	064.8213S	064.8215S																																																																																																																																																	
Diameter (Ø)																																																																																																																																																			
GH	Ø 5.0 mm																																																																																																																																																		
1.5 mm	064.4232S																																																																																																																																																		
Diameter (Ø)																																																																																																																																																			
GH	Ø 5.0 mm																																																																																																																																																		
1.5 mm	064.4233S																																																																																																																																																		

**Auxiliary Instruments**

Torque Control Device		066.1100	
Ratchet		046.119	
Holding Key		046.064	
Implant depth gauge		066.2000	
<b>BL Bone Profiler Kit</b>			
BL Bone Profiler 1		026.0022	
BL Bone Profiler 2		026.0023	
BL Bone Profiler 3		026.0024	
Guiding Pin for BLX Implants		066.0025S	
<b>Abutment removal kit</b>			
RB/WB Abutment-Removal Screw		to push out Abutments, which cannot be removed from the implant due to friction, TAN for BLX Basal Screw, left-hand, length 27 mm, stainless steel	065.0007
Removal Tool		for BLX Basal Screw, left-hand, length 27 mm, stainless steel	065.0008
		for BLX Basal Screw, left-hand, length 21 mm, stainless steel	065.0009
<b>Selection guide for the Straumann® Modular Cassette and instruments</b>			
The Straumann® Modular Cassette Selection Guide helps to define the required Modular Cassette parts and instruments needed for a specific indication. Each chapter includes an order list as well as a set-up overview which shows how to assemble the parts. Art. Nr. 707824			

**Implant-abutment compatibility**

RB/WB and WB		RB Implants Ø 3.5 - 4.5 mm	WB Implants Ø 5.0 - 6.5 mm
«RB/WB» RB/WB Variobase® for Crown		compatible with 	
«WB» WB Variobase® for Crown			

**BLX Healing Abutment selection Quick Guide**


The BLX Healing Abutment selection determines the available options for the final restoration. Therefore, always communicate the implant diameter and the selected healing abutment to the restoring dentist or dental lab in order to facilitate the appropriate final abutment selection.  
Straumann® BLX Implant System – Healing Abutment selection Quick Guide (704658/en).



**International Headquarters**  
Institut Straumann AG  
Peter Merian-Weg 12  
CH-4002 Basel, Switzerland  
Phone +41 (0)61 965 11 11  
Fax +41 (0)61 965 11 01  
www.straumann.com

**RB/WB prosthetics – fit on all BLX implants (∅ 3.5 mm - ∅ 6.5 mm)**

- crowns with slim emergence profiles
- bridges



**WB prosthetics – fit only on WB implants (∅ 5.0 mm - ∅ 6.5 mm)**

- crowns with wide emergence profiles



**Abutment Level**

Implant Analog for Screw-retained Abutments, ∅ 4.6 mm, 0°	Implant Analog for Screw-retained Abutments, ∅ 4.6 mm, 17°/30°	Open-tray Impression for Screw-retained Abutments, ∅ 4.6 mm	Closed-tray Impression for Screw-retained Abutments, ∅ 4.6 mm
025.0050* *for edentulous cases	023.4757	<b>Crown</b> 025.2244 <b>Bridge/Bar</b> 025.0012	<b>Crown</b> 025.2246 <b>Bridge/Bar</b> 025.0014

**Implant Level**

Analog	Open-tray	Closed-tray
065.0021	065.0031 065.0033	065.4310

**Implant Level**

Analog	Open-tray	Closed-tray
065.0022	065.0032 065.0034	065.4810

CARES® Scanbody for Screw-retained Abutments, ∅ 4.6 mm	Repositionable Implant Analog for Screw-retained Abutments, ∅ 4.6 mm	Polishing Aid for Screw-retained Abutments, ∅ 4.6 mm
025.0001	025.0008	025.0005

Straumann® Scanbody	Repositionable Implant Analog	ScanPost for CEREC™ workflow
065.0035	065.0023	065.0038

Repositionable Implant Analog
065.0024

**Temporary Titanium Coping for Screw-retained Abutments, ∅ 4.6 mm**

	∅ 4.6 mm (RB/WB)	
	Crown	024.0023
Replacement screw	Bridge	024.0024
		023.4763

**RB/WB Temporary Abutment, VITA CAD-Temp®**

	GH	1.5 mm	∅ 3.8 mm (RB/WB)	064.4361	064.4371
			∅ 4.5 mm (RB/WB)		

**Temporary Abutments for Crown, AH 10 mm**

	GH	1.5 mm	∅ 3.8 mm (RB/WB)	064.4362	064.4372	–
			∅ 4.5 mm (RB/WB)	064.4363	064.4373	064.4382
			∅ 6.0 mm (RB/WB)	064.4364	064.4374	064.4383
Replacement screw				065.0036		

**Temporary Abutments for Bridge, AH 10 mm**

	GH	1.5 mm	∅ 4.5 mm (RB/WB)	064.4352
				065.0036

**Immediate Temporary Abutment for Crown, AH 5.5 mm**

	GH	1.5 mm	∅ 3.8 mm (RB/WB)	064.43225	064.43325
			∅ 4.5 mm (RB/WB)	064.43235	064.43335
			∅ 5.5 mm (WB)	064.43245	064.43345
Replacement screw				023.0033V2	

**WB Temporary Abutment, VITA CAD-Temp®**

	GH	1.5 mm	∅ 5.5 mm (WB)	064.4390
				065.0036

**Temporary Abutments for Crown, AH 11 mm**

	GH	0.75 mm	∅ 5.5 mm (WB)	064.4391
			1.5 mm	064.4392
				065.0036

**Screw-retained Abutments**

	GH	Angle	∅ 4.6 mm (RB/WB)			
			0.75 mm	0°	17°	30°
			1.5 mm	062.47225	–	–
			2.5 mm	062.47235	–	–
			3.5 mm	062.47245	062.47335	062.47435
			4.5 mm	062.47255	062.47345	062.47445
5.5 mm	–	062.47355	062.47455			
Replacement screw			065.0036			

**RB/WB SRA PLAN Abutment, POM**

	0°	17°	30°	
		025.0073V4	025.0074V4	025.0075V4

**RB/WB Pre-milled Abutment Blanks**

For Medentika Holder, ∅ 11.5 mm	For Medentika Holder, ∅ 15.8 mm	For M-Series, ∅ 12 mm
062.4601	062.4602	062.4603

**Variobase® for Crown, AH 5.5 mm**

	GH	1.5 mm	∅ 3.8 mm (RB/WB)	062.4934	062.4944
			∅ 4.5 mm (RB/WB)	062.4935	062.4945
			2.5 mm	062.4936	062.4946
			3.5 mm	062.4936	062.4946
Replacement screw				065.0036	

**Variobase® for Bridge/Bar Cylindrical, AH 3.5 mm**

	GH	1.5 mm	∅ 4.5 mm (RB/WB)	062.4961
				065.0036

**Variobase® for Crown AS, AH 5.5 mm**

	GH	1.5 mm	∅ 4.5 mm (RB/WB)	062.4972
				065.0037

**Variobase® C**

	GH	1.5 mm	∅ 3.8 mm (RB/WB)	062.4981	062.4982
			∅ 4.5 mm (RB/WB)		
Replacement screw				065.0036	

**Variobase® for Crown, AH 5.5 mm**

	GH	0.75 mm	∅ 5.5 mm (WB)	062.4953
			1.5 mm	062.4954
				065.0036

**Variobase® for Crown AS, AH 5.5 mm**

	GH	1.5 mm	∅ 5.5 mm (WB)	062.4971
				065.0037

**Variobase® C**

	GH	1.5 mm	∅ 5.5 mm (WB)	062.4983
				065.0036

**Burn-out Copings, Plastic**

	∅ 4.6 mm	Crown	023.4748
		Bridge	023.4758
		incl. screw	023.4763

**Variobase® for Bridge/Bar Cylindrical**

	∅ 4.6 mm	Crown	023.0028
		Bridge	023.4754
		incl. screw	023.4763

**Gold Copings**

	∅ 4.6 mm	Crown	023.4753
		Bridge	023.4754
		incl. screw	023.4763

**Occlusal Screw**

	023.4763
--	----------

**CARES® Screw-retained Bars and Bridges**

**CARES® Abutments (RB/WB)**

Individually CAD/CAM manufactured. Please check portal to see options.

**Anatomic Abutments**

	GH	Angle	∅ 6.0 mm (RB/WB)	
			0°	17°
			2.5 mm	062.4103
3.5 mm	062.4104	062.4154		
Replacement screw			065.0036	

**Gold abutments**

	GH	1.5 mm	∅ 3.8 mm (RB/WB)	062.4410	062.4420 (crown)	062.4430 (bridge)
			∅ 4.5 mm (RB/WB)			
Replacement screw				065.0036		

**Novaloc®**

	GH	Angle	∅ 3.8 mm (RB/WB)		
			0°	15°	
			1.5 mm	062.4501	–
			2.5 mm	062.4502	062.4507
			3.5 mm	062.4503	062.4508
			4.5 mm	062.4504	062.4509
5.5 mm	062.4505	062.4510			
6.5 mm	062.4506	062.4511			
7.5 mm	–	062.4512			
Replacement screw			065.0036		

**Processing Package titanium**

	2010.601-NOV
--	--------------

**Processing Package PEEK**

	2010.611-NOV
--	--------------

**CARES® Screw-retained Bars and Bridges**

**WB**

Please check portal to see options.

**CARES® Abutments**

**WB**

Individually CAD/CAM manufactured. Please check portal to see options.

**Gold abutments**

	GH	1.5 mm	∅ 5.5 mm (WB)	062.8410
				065.0036