

STRAUMANN® iEXCEL
EMPOWERING
YOU TO EXCEL



STRAUMANN® iEXCEL SPARK YOUR EXCELLENCE

Straumann® iEXCEL is our ultimate Dental Performance System. Versatility of four implant designs in a single system. Simplicity of one prosthetic connection. Integrated digital workflows and services that help save your time.

* TRUE CLINICAL VERSATILITY

With four implant designs in one single instrument set, you're able to choose the right option to tackle most patient situations, even during surgery. The SLActive® surface excels even for patients who are compromised by smoking, diabetes, or radiation therapy.

* SIMPLICITY YOU'LL LOVE

One TorcFit™ connection for all four implant designs sharing the same connection diameter for all prosthetics. Profit from restorative simplicity with a complete range of Straumann® original stock and customized prosthetics.

* CONNECTED DENTISTRY TO EMPOWER EXCELLENCE

Our system offers integrated digital workflows and services designed to deliver predictable, efficient and high-quality outcomes, from single to full-arch indications.

84%

of clinicians who tried the iEXCEL system would switch to it*



“
The possibility to prepare an osteotomy for both TLX and TLC implants at the same time with the same instrument kit.”

Dr. Simone Janner, Oral Surgeon, Switzerland



*n=50, Jan 2024, Satisfaction survey with clinicians from Europe.

ONE-FOUR-ONE

TRUE CLINICAL VERSATILITY

ONE SET OF INSTRUMENTS

Now you can choose your implant option, even during surgery. A single instrument set works with all four Straumann® iEXCEL implant designs.

FOUR IMPLANT DESIGNS

Address almost every patient situation and excel in immediate treatment protocols with four implant designs in a single system.

BLX AND TLX: EXCEL IN IMMEDIACY


Featuring Dynamic Bone Management, our unique implant design makes immediate protocols achievable and predictable – even in challenging extraction sites, defective sites, or soft bone.

BLC AND TLC: CONFIDENCE TO PERFORM IN ALL INDICATIONS

Our new generation of tapered implants with slim apexes are designed for less invasive treatments, efficient insertion and high primary stability.

ONE SIZE PROSTHETIC CONNECTION

Streamline prosthetic stock thanks to same prosthetic connection across all implant diameters.





Straumann® SLActive® surface demonstrates 100 % implant survival rate, even in patients compromised by smoking, diabetes, or radiation therapy.¹⁻³



Straumann® Roxolid® alloy is stronger and biologically superior to titanium, enabling less invasive and bone preserving treatment options using smaller, stronger implants.⁴⁻⁷





CONFIDENCE IN ALL INDICATIONS

C-DESIGN

- Our new generation of tapered implants allows less invasive treatments. Reduced **3.75 mm diameter** works **in all tooth positions** thanks to the combined use of the Roxolid® material and of the TorcFit™ prosthetic connection.
- With their **slim apex**, the BLC and TLC implants are the option of choice for narrow interdental spaces interventions and for patients with limited anatomy.
- The double-lead thread provides a **faster insertion** while the extended cutting flutes enhance the implant placement by removing the need for tapping.

EXCEL IN IMMEDIACY

X-DESIGN

- Our fully tapered implant design featuring **Dynamic Bone Management**, makes immediate protocols achievable and predictable, even in challenging situations.
- **The extra-short implant BLX 6 mm** is a reliable treatment option in cases where bone height is minimal and can help avoid sinus lift or vertical augmentation.
- The implant's cutting properties **provide direction control** during insertion, helping you achieve optimal positioning.

CHOOSE. MIX. EXCEL

BONE LEVEL

Straumann® Bone Level Implants are designed for a natural look and feel, with Consistent Emergence Profiles™ soft-tissue management providing outstanding esthetic results. Now with Straumann® iEXCEL you can choose from two bone level designs to tackle most patient situations intraoperatively.

BLC and BLX Ø 3.75 mm can be used in all tooth positions to allow bone preservation and highly esthetic outcomes.

BLC FEATURES

1 **PRECISE CORTICAL PREPARATION**

Perfectly matching profile drills for accurate preparation of the crestal bone.

2 **TAPERED DESIGN**

Enables under-preparation of the bone and ensures high implant stability.

3 **EXTENDED CUTTING FLUTES**

Self-tapping properties remove the need for tapping.

4 **SLIM APEX**

For patient's limited anatomy, such as converging root tips, facial undercuts or narrow atrophied ridges.

BLX FEATURES

1 **REDUCED NECK DIAMETER**

Reduces stress on crestal bone and facilitates subcrestal implant placement.

2 **VARIABLE THREAD DESIGN**

Thin and progressive thread design for high implant stability.

3 **BI-DIRECTIONAL CUTTING ELEMENTS**

Enables direction control during insertion, to achieve prosthetically optimal positioning.



4 **DEEP APICAL THREADS**

Large anchoring surface for immediate engagement and improved stability.

	Ø 3.3 mm	Ø 3.5 mm	Ø 3.75 mm	Ø 4.0 mm	Ø 4.5 mm	Ø 5.0 mm	Ø 5.5 mm	Ø 6.5 mm
BLC	✓		✓		✓		✓	✓
BLX		✓	✓	✓	✓	✓	✓	✓
Lengths	8–18 mm	8–18 mm	6–18 mm	6–18 mm	6–18 mm	6–18 mm	6–16 mm	6–16 mm

BLC IMPLANT

BLX IMPLANT



SLActive®
Roxolid®

CHOOSE.MIX.EXCEL

TISSUE LEVEL

The Straumann® iExcel Tissue Level Implants respect the key biological principles of hard and soft tissue healing. They are designed to significantly reduce the risk of inflammation and bone resorption as the implant-abutment interface is moved away from the bone.⁸ With the prosthetic platform above bone level they help to achieve better visibility and accessibility, even in the posterior region.

Tissue level implants are a good choice for patients with a history of periodontal conditions.

TLC FEATURES

1 PRECISE CORTICAL PREPARATION

Perfectly matching profile drills for accurate preparation of the crestal bone.

2 TAPERED DESIGN

Enables under-preparation of the bone and ensures high implant stability.

3 EXTENDED CUTTING FLUTES

Self tapping properties remove the need for tapping.

4 SLIM APEX

For patient's limited anatomy, such as converging root tips, facial undercuts or narrow atrophied ridges.

TLX FEATURES

1 REDUCED NECK DIAMETER

Reduces stress on crestal bone and facilitates subcrestal implant placement.

2 VARIABLE THREAD DESIGN

Thin and progressive thread design for high implant stability.

3 BI-DIRECTIONAL CUTTING ELEMENTS

Enables direction control during insertion to achieve prosthetically optimal positioning.



4 DEEP APICAL THREADS

Large anchoring surface for immediate engagement and improved stability.

	Ø 3.3 mm	Ø 3.5 mm	Ø 3.75 mm	Ø 4.0 mm	Ø 4.5 mm	Ø 5.0 mm	Ø 5.5 mm	Ø 6.5 mm
TLC	✓		✓		✓		✓	✓
TLX			✓		✓		✓	✓
Lengths	8–18 mm	–	6–18 mm	–	6–18 mm	–	6–16 mm (TLC) 6–12 mm (TLX)	6–16 mm (TLC) 6–12 mm (TLX)

TLC IMPLANT

TLX IMPLANT



SLActive®

Roxolid®

MODULAR CASSETTE AND VELODRILL™ SYSTEM

ONE INSTRUMENT SET

* STRAUMANN® MODULAR CASSETTE

The Straumann® Modular Cassette is a compact instrument set which can be assembled and combined in different ways, allowing you to customize your kit with the tools you need.

- Module A: for universal and prosthetic instruments
- Module B: for workflow and drills – one insert that works with all iExcel implants
- Module C: add-on for guided surgery (handles) or a simple box for auxiliaries

* STRAUMANN® VELODRILL™ FLEXIBLE. SIMPLE. COOL.

The VeloDrill™ system enables you to use a single drill set that works with all Straumann® iEXCEL implant designs.

FLEXIBLE


Shared instrumentation for freehand and guided surgery with both X and C designs.

SIMPLE

Drills are color-coded for easy identification of diameters.


COOL

Low-temperature drilling. Avoids overheating surrounding structures.



NO LIMITS. JUST CHOICE.

SIMPLICITY YOU'LL LOVE





“
After 5 years of using BLX, I can say
that the TorcFit™ connection is awesome –
crestal bone stability is outstanding. ”

Dr. Algirdas Puisys, Periodontologist, Lithuania

* TORCFIT™ CONNECTION

STRENGTH OF TORX INDEX

- Allows transfer of high torques
- Six positions for flexible implant-abutment alignment

PRECISE FIT OF 7° CONICAL CONNECTION

- Precise implant-abutment fit
- Abutment screw only engages in implant when the abutment is correctly seated
- Correctly seated abutments are secured with friction fit
- High mechanical stability and stress distribution

ONE SIZE PROSTHETIC CONNECTION ACROSS ALL IMPLANTS

- One single driver for all implants
- Lower prosthetic inventory
- Simpler staff training

EXCELLENCE THROUGH COMPLETE PROSTHETIC PORTFOLIO

ALL INDICATIONS AND DIGITAL WORKFLOWS.

Consistent Emergence Profiles™
throughout the full portfolio,
for easy soft-tissue conditioning
across all indications.





STRAUMANN® VARIOBASE®

Our Variobase® portfolio provides efficient production without compromising restorative flexibility. Available with angled screw-channel options and in various heights for single-tooth to edentulous indications. The specially designed rotational locks offer an increased bonding surface for maximum reliability.



PRE-MILLED ABUTMENT BLANKS

Pre-milled abutment blanks enable in-house manufacturing of one-piece custom abutments. With a pre-machined original Straumann® connection, precision and reliability are assured.



CUSTOMIZED PROSTHETICS

Customized Prosthetics ensure ultimate restorative outcomes, with the abutment emergence profile and shape specifically designed for your patient. The screw-retained bar and bridge offer angle correction to solve the most demanding cases.



SCREW RETAINED ABUTMENTS (SRA) FOR BRIDGES AND FULL ARCH

The smooth, cleanable and undercontoured design reduces the need for bone profiling in surgery and is designed to maximize space for soft tissue.



STRAUMANN® NOVALOC® RETENTIVE SYSTEM

Low friction between the abutment and matrix ensures excellent durability for removable indications. The 15° angled abutments reduce unilateral stress and wear.



EXCEL WITH THE STRAUMANN® DIGITAL WORKFLOW FROM SINGLE UNIT TO FULL-ARCH

The Straumann® digital workflow empowers your growth with a comprehensive range of straightforward and open solutions, streamlining and upgrading your workflow for the final benefit of your patients. You can confidently rely on a single provider for all your implant procedures.

SCANNING

Straumann SIRIOS™



PATIENT ENGAGEMENT


Smilecloud



smilecloud

TREATMENT PLANNING

coDiagnostiX®



Smile in a Box®

smilecloud

INTERDISCIPLINARY COLLABORATION

STRAUMANN SIRIOS™*


Go beyond image capturing and open the door to digital dentistry with Straumann SIRIOS™; with ultimate accuracy and a revolutionary scanning performance, it is specifically designed to unlock streamlined digital workflows.

SMILECLOUD

Captivate your patients through natural biometric smile design and effectively engage your network with the potential of interdisciplinary workflows. Start today with Smilecloud!

CODIAGNOSTIX®

Unlock the potential of treatment planning with coDiagnostiX® including AI-powered data preparation, visualization and automated tooth extraction designed to efficiently access immediate treatments.



Discover the Straumann® solutions for digital workflow.



FLEXIBILITY

From single unit to edentulous cases, with our solutions you have the right tools to streamline your workflow.



PREDICTABILITY

Boost your confidence with digital treatment planning and guided procedures designed to achieve predictable results.



EFFICIENCY

Implement streamlined workflows that are designed to promote speed and quality in your treatment.



EXCELLENCE POWERED BY AI

Our solutions are infused with AI to empower you to achieve your goals with efficiency and speed.

SURGERY

SECONDARY CANNING

CUSTOMIZED PROSTHETIC DESIGN & MANUFACTURE

Straumann® iEXCEL Performance System



Straumann SIRIOS™



RevEx™



Straumann® UNIQ



SMILE IN A BOX®

Outsource planning, design, and/or manufacturing and access immediate treatment protocols without requiring additional hardware and software investments.

STRAUMANN® REVEX™ SCANBODY

Capture full-arch implant positions to enable a full-arch workflow from temporary to final restoration. RevEX™ provides greater patient comfort, with fewer visits.⁹⁻¹¹

STRAUMANN® UNIQ*

Use a premium service to outsource the planning, design, and manufacturing of your custom implant prosthetics on demand, based on your specific needs.

*The availability varies by region; contact local team for details.

STRAUMANN® PRO ARCH*

EXCEL IN FULL-ARCH

Powered by iEXCEL, our ultimate dental performance system, Straumann® Pro Arch now delivers more performance and versatility with four implant designs in a single instrument set. Evidence-based technologies such as Roxolid® and SLActive®, and integrated digital workflows provide clinicians with predictability, and peace of mind in challenging clinical situations while bringing efficiency to your practice.

* PATIENT-CENTRIC VERSATILITY

Bone and Tissue Level implants

→ Length: from 6 mm to 18 mm

→ Diameter: from Ø3.3 mm to Ø6.5 mm



SLActive®



Roxolid®



* SIMPLICITY FOR CHALLENGING CASES

→ One connection

→ One instrument set

→ Short surgical protocols

* DIGITAL FULL-ARCH

→ Straumann® Smile in a Box®

→ Straumann® RevEX™



1 NARROW IMPLANTS

Ø3.3 mm BLC and Ø3.5 mm BLX implants with Roxolid® to preserve the bone and avoid bone grafting.

2 SHORT IMPLANTS

BLX 6 mm can help avoid sinus lift or vertical augmentation in cases where bone height is minimal.

3 SMALL ALL INDICATION IMPLANT

Ø3.75 mm for all tooth positions to reduce invasiveness and preserve bone.

4 ROUNDED ANGULATED ABUTMENT

Designed to maximize space for soft tissue.


5 STRAUMANN® UNIQ*

Comprehensive bar portfolio.

*The availability varies by region; contact local team for details.

STRAUMANN® REGENERATIVE SOLUTIONS MASTER ANY CHALLENGE

To ensure clinical success, dental implants require stable and healthy hard and soft tissue. Straumann® Regenerative Solutions regenerate the hard and soft tissue around teeth or dental implants, to maintain long-lasting functionality and esthetics. Our solutions are built on three pillars along the patient's dental journey, and cover the entire range of cases, from simple to complex.





TOOTH PRESERVATION

Natural teeth are unique and, as dental professionals, we strive to preserve them as long as possible. Untreated inflammation can advance to periodontitis and potentially result in tooth loss.

Straumann® Emdogain® is a central part of our solution to support the effective regeneration of the periodontium and accelerate the healing to help maintain natural teeth.



IMPLANT SITE MANAGEMENT

The right amount of hard and soft tissue is key for successful dental implant placement. Our solution includes a comprehensive biomaterials portfolio to support tissue regeneration for optimal implant placement conditions.



IMPLANT PRESERVATION

Healthy peri-implant tissues are vital to preserve the placed implant. Untreated inflammation can advance to peri-implantitis and potentially result in implant loss. Straumann® GalvoSurge® is a central part of our solution to support the effective regeneration of hard and soft tissue to help maintain dental implants.

GALVOSURGE® DENTAL IMPLANT CLEANING SYSTEM

- **Speed:** effectively removes biofilm from the implant surface in just 2 minutes
- **Reliability:** provides optimal conditions for bone regeneration
- **Safety:** offers a minimally invasive treatment process to successfully treat peri-implantitis without implant removal



CLINICAL CASES

Lower molar tooth replacement in healed site


Freehand surgery for a lower molar tooth replacement in a healed site using an apically tapered Straumann® BLC implant Ø4.5 mm × 12 mm SLActive®, Roxolid®, and following a conventional loading protocol.



Dr. Algirdas Puišys



Intraoral situation of left lower quadrant




CBCT analysis reveals available bone before implantation



Implant bed ready to receive the implant



Straumann® BLC implant Ø4.5 mm × 12 mm, SLActive®, Roxolid® in final position




Post insertion radiograph



12 weeks of osseointegration and soft tissue healing



Definitive screw-retained implant crown after 12 weeks



Radiograph of final crown on Variobase® abutment screw-retained in final position

Immediate implant placement and provisional restoration of upper premolar


Freehand surgery for immediate implant placement of an upper premolar, using an apically tapered Straumann® BLC implant Ø4.5 mm × 12 mm, SLActive®, Roxolid®, and immediate provisional restoration, followed by final restoration using a digital workflow.



Dr. Algirdas Puišys



Intraoral situation of tooth #24



CBCT analysis



Occlusal view of the extraction socket



Insertion of Straumann® BLC implant Ø4.5 mm × 12 mm, SLActive®, Roxolid®



Implant in final position prior to gap management with bone grafting



Immediate provisional crown



12 weeks of osseointegration and soft tissue healing




Definitive screw-retained implant crown after 12 weeks

Lower molar tooth replacement in healed site

Utilizing digital planning for the replacement of a lower molar tooth in a healed site to establish a prosthetically-driven implant position. Conducting free-hand surgery for the placement of an apically tapered Straumann® TLC implant, Ø4.5mm × 10 mm SLActive®, Roxolid®. Following a healing period of 14 weeks, the digitally fabricated screw-retained crown is delivered.




Dr. Alessandro Perucchi



Dr. Cristian Scognamiglio



Intraoral situation of right lower quadrant



Digital planning to establish a prosthetically-driven implant position



Full-thickness flap was raised to allow access to the bone and evaluate the local anatomy. Alignment pin placed to check the 3D position of the osteotomy and preparation depth.



Implant in final position



Radiographical control after insertion of dental implants.



Scanbody for digital impression



Definitive screw-retained crown in position after 14 weeks




Radiograph of final crown on Variobase® abutment screw-retained in final position

Maxillary overdenture with immediate loading

Edentulous maxilla with limited bone height in the posterior region. Flapless guided surgery, placement of two Straumann® BLC implants Ø3.75 × 10 mm (at positions 12 and 22) and two Straumann® BLC implants Ø4.5 × 10 mm (at positions 14 and 24). Immediate fixed temporization was achieved using the patient's former denture, which remained in function for three months. Subsequently, considering the patient's presentation with a thin upper lip, high smile line, and vertical bone crest atrophy, a final metal-reinforced overdenture fixed on a Dolder bar was selected.




Dr. Jiri Krug



Occlusal view before surgery



Fully guided drilling



Implant placement, BLC Ø3.75 × 10 mm SLActive®, Roxolid®



Provisional bridge with the titanium copings made from a reinforced full denture



Screw-Retained Abutments after 8 weeks



Position of the CAD/CAM titanium Dolder bar (Createch) within the outline of the overdenture



Panoramic X-ray done for Sheffield test approval



Post-op smile at 4-month follow up

DRIVEN TO EXCEL

FOR YOU AND YOUR PATIENTS

A comprehensive array of products, solutions, and services made Straumann® a leader in the global market for implant, restorative, and regenerative dentistry. With a focus on digital, seamless workflows, Straumann® provides an end-to-end, customer-centric treatment experience. An integrated ecosystem enhances efficiency, convenience, and collaboration, supporting dental professionals with secure data storage and sharing, and leveraging evidence-based and data-driven processes.

Available in more than 100 countries, the Straumann® brand is associated with Swiss quality, precision, and innovation, fostering confidence among dental professionals, supported by the world's largest global scientific network and a partnership with the International Team for Implantology (ITI).



Straumann® is the flagship brand of the world's largest implant manufacturer group – the Straumann Group. Renowned as one of the most iconic brands in dentistry, we have a long-standing legacy of pioneering significant advancements in the field. We are at the forefront of expanding the scope of dental care by relentlessly and consistently providing excellence and transformational innovation to a global community of forward-thinking dental professionals.



STRAUMANN® iEXCEL – EMPOWERING YOU TO EXCEL



Discover
Straumann® iEXCEL

- Straumann® iEXCEL is our ultimate Dental Performance System
- Versatility of four implant designs in a single system
- Simplicity of one prosthetic connection
- Integrated digital workflows and services that help save your time

REFERENCES

1 Alsahhaf A et al. Survival of Titanium-Zirconium and Titanium Dental Implants in Cigarette-smokers and Never-smokers: A 5-Year Follow-up. Chin J Dent Res. 2019;22(4):265-272 2 Cabrera-Domínguez JJ, Castellanos-Cosano L, Torres-Lagares D, Pérez-Fierro M, Machuca-Portillo G. Clinical performance of titanium-zirconium implants with a hydrophilic surface in patients with controlled type 2 diabetes mellitus: 2-year results from a prospective case-control clinical study. Clin Oral Investig. 2020 Jul;24(7):2477-2486 3 Nelson, K., Stricker, A., Raguse, J.-D. and Nahles, S. (2016), Rehabilitation of irradiated patients with chemically modified and conventional SLA implants: a clinical clarification. J Oral Rehabil, 43: 871 -872 4 Bernhard N et al. Forum Implantologicum 2009; 5(1): 30-39 5 Hotchkiss KM et al. Clin Oral Implants Res. 2017 Apr;28(4):414-423 6 Al-Nawas B. et al. J Oral Implantol. 2015 Aug;41(4):e118-25 7 Chiapasco M. et al. Clin Oral Implants Res. 2012 Oct;23(10):1136-41 8 Derks J, Schaller D, Hakansson J, Wennstrom JL, Tomasi C, Berglundh T. Effectiveness of Implant Therapy Analyzed in a Swedish Population: Prevalence of Peri-implantitis. J Dent Res. 2016 Jan; 95(1):43-9. 9 Panos Paspaspyridakos, Armand Bedrossian , Yukio Kudara , Panagiotis Ntovas , Abdullah Bokhary , Konstantinos Chochlidakis . Reverse scan body: A complete digital workflow for prosthesis prototype fabrication. Journal of Prosthodontics, 2022 10 E. Armand Bedrossian DDS, MSD, FACP, Panos Paspaspyridakos DDS, MS, PhD, Edmond Bedrossian DDS, FACOMS, FAO, FITI “The Reverse Scan Body protocol”: Completing the Digital Workflow, Compendium June 2023 (will be published in August) 11 Vandeweghe S, Vervack V, Dierens M, De Bruyn H. Accuracy of digital impressions of multiple dental implants: an in vitro study. Clin Oral Implants Res 2016;28:648-53.

International Headquarters

Institut Straumann AG
Peter Merian-Weg 12
CH-4002 Basel, Switzerland
Phone +41 (0)61 965 11 11
Fax +41 (0)61 965 11 01
www.straumann.com

© Institut Straumann AG, 2024. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

 **straumann**