

DIGITAL MATERIALS OVERVIEW

MATERIALS ¹

	MILLING				SERVICE / SOLUTION			APPLICATION ¹										
	CARES® C Series	CARES® M series	CARES® D series	Centralized milling	CARES® X-Stream™	CARES® Scan and Shape	Connectivity	Tooth restoration			Others	Direct connection to implants						
								Inlays ² / Onlays ² / Veneers	Copings / Crowns	Bridges / Bars (maximum units)	Splints	CARES® Customized Abutments and Pre-Milled Blanks	etkon® iDent Customized Abutments and Pre-Milled Blanks	CARES® Screw-Retained Bridges and Bars				
Glass-ceramics	Translucency																	
nIce®				●	●		●	●	●									
IPS e.max® CAD	● ³	● ³		●	●		●	●	●									
VITA® Suprinity				●	●													
VITA® Mark II				●	●			●	●									
VITA® TriLuxe / TriLuxe Forte				●	●			●	●									
Hybrids	Translucency																	
3M™ ESPE™ Lava™ Ultimate		●						●	●									
VITA® Enamic	●	●																
Ceramics	Translucency																	
zerion® UTML							●	●	●	●	●	3						
zolid® SHT ML		●										3						
zolid® SHT / SHT preshade		●										3						
zerion® ML							●	●	●	●	●	16						
3M™ ESPE™ Lava™ Plus Zirconia							●	●	●	●	●	16						
zolid® HT / HT preshade		●										16						
zerion® LT							●	●	●	●	●	16						
zi		●										16						
zerion® GI							●	●	●	●	●	16						
Zirconium Dioxide for CARES®							●	●	●	●	●	16		● ⁴				
Metal																		
coron® Cobalt-Chromium							●	●	●	●	●	16		●				●
Cobalt-Chromium for etkon® iDent	●	●																●
Sintron®		●																
TAN for CARES®	●	●						●	●	●	●			●				
ticon® Titanium Grade 4								●	●	●	●							●
Titanium for CARES®							●	●	●	●	●			● ⁴				
Titanium for etkon® iDent	●	●						●	●	●	●							●
Polymer																		
polycon® ae							●	●	●	●	●	16						
VITA® CAD-Temp	●	●																
Temp / Temp multilayer		●																
PEEK		●																
Splintec		●												●				
polycon® cast								●	●	●	●	16						
PMMA (burn out polymer)		●																
Wax		●																

¹ application and material availability might differ from country to country.

² veneerable inlay/onlay for 3M™ ESPE™ Lava™ Plus Zirconia, zerion®, ticon®, coron®, polycon® ae and polycon® cast materials.

³ material listed and mill-able

⁴ CARES® Abutments available, no Pre-Milled Abutments Blanks available

MO: Medium Opacity LT: Low Translucency HT: High Translucency GI: Gingiva shaded

VITA® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Germany.

3M™, ESPE™ and Lava™ are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada.

JUVORA™ is a trademark of JUVORA Ltd.

IPS e.max® is a registered trademark of Ivoclar Vivadent AG, Liechtenstein.