



Components overview  
2021 - 2022



AXIOM® BONE LEVEL SURGICAL COMPONENTS		AXIOM® BONE LEVEL SURGICAL COMPONENTS								AXIOM® BONE LEVEL SURGICAL COMPONENTS																			
IMPLANT RANGE		AXIOM® BL REG				AXIOM® BL PX				HEALING SCREW	HEALING SCREWS																		
STERILE Cover screw included OPIM100		IMPLANTS								STERILE Manual	AXIOM® BL REG & PX																		
											Ø 3.4 mm				Ø 4.0 mm				Ø 5.0 mm				Ø 6.0 mm						
											Coronal height 2 mm		1 mm		2 mm		1 mm		2 mm		1 mm		2 mm		1 mm				
mm											Gingival height																		
6.5											0.75 mm																		
8											1.5 mm																		
10											2.5 mm																		
12											3.5 mm																		
14											4.5 mm																		
16																													
18																													
Round bur* (INF20)											Wrench																		
Pointer drill (OPPO15150)											Manual hexagonal surgical wrench																		
Lindemann bur* (OPR20)											OPCS100																		
Initial drill Gauge																													
Ø: mm																													
Stepped drill Gauge																													
Ø: mm																													
	Ø 1.5																												
	Ø 2.0																												
	Ø 2.0																												
	Ø 2.4																												
	Ø 3.0																												
	Ø 3.6																												
	Ø 4.2																												
	Ø 4.8																												
	Ø mm																												
	Implant mandrel / Wrench	SHORT		MEDIUM		LONG		SURGICAL DYNAMOMETRIC WRENCH																					
		23 mm	21 mm	27 mm	27 mm	32 mm	31 mm																						
		OPMV180	OPCV060	OPMV215	OPCV110	OPMV250	OPCV160		INCCDC																				

\*Not included in INMODOPS3 surgery kit

**Axiom® BL countersinks for optimized bone crest preparation**

Countersinks with pin (included)	Ø 4.5 mm	OPFF45
	Ø 5.3 mm	OPFF53
	Ø 6.6 mm	OPFF66



AXIOM® BONE LEVEL PROSTHETIC COMPONENTS			
IMPRESSION			
		IMPLANT LEVEL	
		CONVENTIONAL	DIGITAL
		Short	Long
Pop-in transfer	Manual	OPPI100S (x1) OPPI100S-4 (x4)	OPPI100 (x1) OPPI100-4 (x4)
Pick-up transfer	Manual	OPPU100 (x1) OPPU100-4 (x4)	OPPU100L (x1) OPPU100L-4 (x4)
Impression ring for pick up transfers		Ø 4.0 mm	Ø 5.0 mm
		OPPROFIL410	OPPROFIL510
		OPPROFIL420	OPPROFIL520
		OPPROFIL430	OPPROFIL530
Analog	Manual	OPPIA100 (x1) OPPIA100-4 (x4)	152-27-PA (x1) 152-27-PA-4 (x4)
			PA-TOOL-01

**AXIOM® BONE LEVEL PROSTHETIC COMPONENTS**

ABUTMENT						
SINGLE & MULTIPLE CEMENT RETAINED RESTORATIONS						
AESTHETIC ABUTMENTS	Angulation	Gingival height	Ø 3.4 mm	Ø 4.0 mm	Ø 5.0 mm	Ø 6.0 mm
 M1.6 Black screw included OPTS160 Extractible with INEXPS/L	0°	0.75 mm	-	OPAT400	OPAT500	-
		1.5 mm	-	OPAT410	OPAT510	OPAT610
		2.5 mm	-	OPAT420	OPAT520	OPAT620
		3.5 mm	-	OPAT430	OPAT530	OPAT630
	7°	4.5 mm	-	OPAT440	OPAT540	OPAT640
		0.75 mm	-	OPAT40-7	OPAT50-7	-
		1.5 mm	OPAT31-7	OPAT41-7	OPAT51-7	-
		2.5 mm	OPAT32-7	OPAT42-7	OPAT52-7	-
	15°	3.5 mm	OPAT33-7	OPAT43-7	OPAT53-7	-
		4.5 mm	OPAT34-7	OPAT44-7	OPAT54-7	-
0.75 mm		-	OPAT401	OPAT501	-	
1.5 mm		OPAT311	OPAT411	OPAT511	OPAT611	
23°	2.5 mm	OPAT321	OPAT421	OPAT521	OPAT621	
	3.5 mm	OPAT331	OPAT431	OPAT531	OPAT631	
	4.5 mm	OPAT341	OPAT441	OPAT541	OPAT641	
	0.75 mm	-	OPAT402	OPAT502	-	

**AXIOM® BONE LEVEL PROSTHETIC COMPONENTS**

**Axiom® Multi-Level® Prosthetic Kit**

INMODOPP3

**Torq Control® Universal torque wrench**

15501

**AXIOM® BONE LEVEL PROSTHETIC COMPONENTS**

TI-BASE AND PREMILLED FOR LABSIDE CAD/CAM		
SINGLE CEMENTED OR SCREW RETAINED RESTORATIONS		
TI-base	FLEXIBASE	CEREC
 M1.6 Black screw included OPTS160 Hexagonal wrench Extractible with INEXPS/L	<b>Coronal height: 3.5 mm</b>	<b>Size L block</b>
	<b>Gingival height</b> Ø 4.0 mm Ø 5.0 mm	<b>Gingival height</b> Ø 5.0 mm
	1.5 mm OPFLEX413 OPFLEX513	1.5 mm OPBASE-C501-L
	2.5 mm OPFLEX423 OPFLEX523	2.5 mm OPBASE-C502-L
	3.5 mm OPFLEX433 OPFLEX533	3.5 mm OPBASE-C503-L
Castable coping	Ø 4.0 mm OPFLEXC403	Ø 5.0 mm OPFLEXC503
PRE-MILLED	PREFACE	TI-FORMS*
	<b>Medentika® holder</b> Ø 11.5 mm Ø 16.0 mm	<b>Amann Girsch holder</b> Ø 11.5 mm
Definitive titanium M1.6 included	OPPM-M-TI-11	OPPM-M-TI-16
OPTS161		OPPM-C-TI-11
Hexagonal wrench		
Extractible with INEXPS/L		

**AXIOM® BONE LEVEL PROSTHETIC COMPONENTS**

TEMPORIZATION	
IMPLANT LEVEL	
Gingival height	Implant Level
0.75 mm	Ø 3.4 mm, Ø 4.0 mm, Ø 5.0 mm, Ø 6.0 mm
1.5 mm	OPTP310, OPTP410, OPTP510, OPTP610
2.5 mm	OPTP320, OPTP420, OPTP520, OPTP620
3.5 mm	OPTP330, OPTP430, OPTP530, OPTP630
4.5 mm	OPTP340, OPTP440, OPTP540, OPTP640

**AXIOM® BONE LEVEL PROSTHETIC COMPONENTS**

INSTRUMENTS & SCREWS	
Instrument	SKU
Short hexagonal wrench 14 mm	INCHCV
Long hexagonal wrench 21.5 mm	INCHELV
XL hexagonal wrench 29 mm	INCHEXLV
Short hexagonal mandrel 20 mm	INMHECV
Long hexagonal mandrel 27 mm	INMHVELV
XL hexagonal mandrel 35 mm	INMHXLV
Torq Control® Dynamometrical wrench	15501
Dynamometrical ratchet wrench	INCCD
Short disconnecter and gripper for Axiom® BL abt 32.4 mm	INEXPS
Long disconnecter and gripper for Axiom® BL abt 35.4 mm	INEXPL
Black M1.6 prosthetic screw	OPTS160
Titanium M1.6 prosthetic screw	OPTS161
M1.6 laboratory screw	OPTS162
Screw driver for SA	SATOOL-01
Wrench for analog for printed model	PA-TOOL-01

**AXIOM® BONE LEVEL PROSTHETIC COMPONENTS**

ABUTMENT					
REMOVABLE OVERDENTURE					
NOVALOC®		LOCATOR®		BALL ABUTMENTS	
<b>1.5 mm</b>	OPNOVA015	<b>1.5 mm</b>	O1697	<b>1.5 mm</b>	AXDA015
<b>2.5 mm</b>	OPNOVA025	<b>2.5 mm</b>	O1698	<b>2.5 mm</b>	AXDA025
<b>3.5 mm</b>	OPNOVA035	<b>3.5 mm</b>	O1699	<b>3.5 mm</b>	AXDA035
<b>4.5 mm</b>	OPNOVA045	<b>4.5 mm</b>	O1016	<b>4.5 mm</b>	AXDA045
<b>5.5 mm</b>	OPNOVA055	<b>5.5 mm</b>	O1017	-	-

**STANDARD ABUTMENTS**

STANDARD ABUTMENTS								
SINGLE & MULTIPLE CEMENT RETAINED RESTORATIONS								
STANDARD ABUTMENTS	Angulation	Gingival height	Coronal Height	Ø 3.4 mm	Ø 4.0 mm	Ø 5.0 mm		
 M1.6 Black screw included OPTS160 Extractible with INEXPS/L	0°	1.5 mm	4 mm	OPST314	OPST414	OPST514		
		2.5 mm	6 mm	OPST316	OPST416	OPST516		
		3.5 mm	4 mm	OPST324	OPST424	OPST524		
		4.5 mm	6 mm	OPST326	OPST426	OPST526		
	15°	1.5 mm	6 mm	-	OPST416-15	OPST516-15	-	
		2.5 mm	6 mm	-	OPST426-15	OPST526-15	-	
		23°	3.5 mm	6 mm	-	OPST436-15	OPST536-15	-
			4.5 mm	6 mm	-	OPST446-15	OPST546-15	-
	Non indexed		1.5 mm	6 mm	-	OPST416-23	OPST516-23	
	Non indexed		2.5 mm	6 mm	-	OPST426-23	OPST526-23	
Non indexed		3.5 mm	6 mm	-	OPST436-23	OPST536-23		

**IMPRESSION**

ABUTMENT LEVEL				
Standard abutment transfer	OPTT100 (x1) / OPTT100-5 (x5)			
Standard abutment Analog	Ø 3.4 mm	Ø 4.0 mm	Ø 5.0 mm	Ø 6.0 mm
Coronal height 4 mm				
Coronal height 6 mm				

**ACCESSORIES & FINAL COPING**

SINGLE & MULTIPLE CEMENT RETAINED RESTORATIONS				
Accessories	Ø 3.4 mm	Ø 4.0 mm	Ø 5.0 mm	Ø 6.0 mm
PEEK protective cap for standard abutments (temporary cement)	OPPC304	OPPC404	OPPC504	OPPC604
PMMA castable coping for standard abutments	OPCA304	OPCA404	OPCA504	OPCA604
Indexed	OPCR304	OPCR404	OPCR504	OPCR604
Non indexed	OPCR306	OPCR406	OPCR506	OPCR606

**MULTIPLE SCREW RETAINED RESTORATIONS**

MULTI-UNIT ABUTMENTS							
SINGLE & MULTIPLE CEMENT RETAINED RESTORATIONS							
Multiple Screw Retained Restorations	Ø 4.8 mm / R			Ø 4.0 mm / N			
 MU M1.6 Black screw included with angulated abutments OPMU160	Gingival height	straight	18° angulated	30° angulated	straight		
		Indexed		Non indexed			
		0.75 mm	OPMU0-0	OPMU18-1-IN	OPMU18-1	OPMU30-1-IN	OPMU30-1
		1.5 mm	OPMU0-1	OPMU18-2-IN	OPMU18-2	OPMU30-2-IN	OPMU30-2
		2.5 mm	OPMU0-2	OPMU18-3-IN	OPMU18-3	OPMU30-3-IN	OPMU30-3
3.5 mm	OPMU0-3	-	-	-	-		
4.5 mm	OPMU0-4	-	-	-	-		

**ABUTMENT LEVEL**

MULTI-UNIT			
Multi-Unit	Conventional	Digital	Laboratory scan adapter
Pop-in transfer	MUT200 (x1), MUT200-4 (x4)	MUNT200 (x1), MUNT200-4 (x4)	151-03-SA0
Pick-up transfer	MUT100 (x1), MUT100-4 (x4)	MUN100 (x1), MUN100-4 (x4)	151-03-DT-MU, 151-04-DT-MUN, 151-03-PA, 151-04-PA
Abutment analog	MUA100 (x1), MUA100-4 (x4)	MUNA100 (x1), MUNA100-4 (x4)	151-03-PA, 151-04-PA

**MULTIPLE SCREW RETAINED RESTORATIONS**

MULTI-UNIT	
Multiple Screw Retained Restorations	Multi-Unit
Protective cap	MUCAP (x1), MUCAP-4 (x4)
Castable coping	MUC300, MUNC300
CoCr castable coping	MUC400, MUNC400
Straight pacific ring	MUNPAC100
Analog pacific	MUNPAC110
Castable pacific	MUNPAC120

**MULTI-UNIT**

MULTI-UNIT	
Multi-Unit	Multi-Unit
Flexibase	MUFLEX, MUNFLEX
Castable coping	MUFLEXC, MUNFLEXC

**MULTI-UNIT**

MULTI-UNIT	
Multi-Unit	Multi-Unit
Titium provisional coping	MUCI00, MUNCI00
Castable coping	MUC200, MUNC200

**INSTRUMENTS & SCREWS**

INSTRUMENTS & SCREWS	
Instrument	SKU
Wrench for straight MU abt 24 mm	MUW100
Mandrel for straight MU abt 5 16 mm / L 25.5 mm	MUM100, MUM100L
Short prehensile tool for MU	MUWS
MU M1.6 Black screw	OPMU160
MU titanium M1.6 laboratory screw	OPMU161, OPMU161-4
MU M1.4 Black screw	MU140
MU M1.4 definitive screw for zirconia prosthesis	MU140Z
MU titanium M1.4 laboratory screw (blue)	MU141
Protective analog Ø 4.8 mm / R	MUA200 (x7), MUA200-4 (x4)
Protective analog Ø 4.0 mm / N	MUNA200
MU laboratory screw	SHORT MUTI01, LONG MUTI02

**INSTRUMENTS**

INSTRUMENTS	
Instrument	SKU
Equipment box with 3 instruments	2010.101-STM
Processing insert removal instrument	2010.731-STM
Retention insert instrument	2010.741-STM
Matrix housing extraction instrument	2010.751-STM