



Components overview



AXIOM[®] TISSUE LEVEL SURGICAL COMPONENTS

AXIOM[®] TISSUE LEVEL SURGICAL COMPONENTS

AXIOM[®] TISSUE LEVEL SURGICAL COMPONENTS

AXIOM[®] TISSUE LEVEL SURGICAL COMPONENTS

IMPLANT RANGE	IMPLANTS												HEALING SCREWS							
	AXIOM [®] TL REG												AXIOM [®] TL PX							
	R PLATFORM 4.8 mm				N PLATFORM 4.0 mm				R PLATFORM 4.8 mm				N PLATFORM 4.0 mm							
STERILE Cover screw included TSS-RO TSS-NO													STERILE Manual							
Length Neck height													Conical				Straight			
1.5 mm													2 mm				2 mm			
2.5 mm													4 mm				4 mm			
3.5 mm													2 mm				2 mm			
5 mm																				
6.5 mm																				
8 mm																				
10 mm																				
12 mm																				
14 mm																				
16 mm																				
Round bur* (INFB20) Painter drill (OPPO15150) Lindemann bur* (OPR20)																				
Initial drill																				
Gauge																				
Stepped drill																				
Gauge																				
Implant tap (optional)																				
Implant mandrel / wrench																				

* Not included in INMODOP53 surgery kit

**Axiom[®] Multi-Level[®]
Surgical kit**

INMODOP53



Manual surgical wrench



OPCS100

AXIOM® TISSUE LEVEL PROSTHETIC COMPONENTS

AXIOM® TISSUE LEVEL PROSTHETIC COMPONENTS

AXIOM® TISSUE LEVEL PROSTHETIC COMPONENTS

AXIOM® TISSUE LEVEL PROSTHETIC COMPONENTS

AXIOM® TISSUE LEVEL PROSTHETIC COMPONENTS

Axiom® Multi-Level® Prosthetic Kit
INMODOPP3

	IMPRESSION TAKING				Wrench for analog for printed model PA-TOOL-01
	SINGLE RESTORATION				
	Conventional impression		Digital impression		
Pop-in transfer Screw included Hexagonal wrench	R PLATFORM 4.8 mm Short	R PLATFORM 4.8 mm Long	N PLATFORM 4.0 mm Short	N PLATFORM 4.0 mm Long	Digital transfer
	TT200-R	TT200L-R	TT200-N	TT200L-N	
Pick-up transfer Indexed Long and short screws included Hexagonal wrench	TT300-R		TT300-N		156-02-DT
	TA100-R (x1) TA100-R-4 (x4)		TA100-N (x1) TA100-N-4 (x4)		156-01-PA (x1) 156-02-PA-4 (x4)
Abutment analog Hexagonal wrench	TA100-R (x1) TA100-R-4 (x4)		TA100-N (x1) TA100-N-4 (x4)		156-02-PA (x1) 156-01-PA (x1)

	MULTIPLE SCREW RETAINED PROSTHESIS				Wrench for analog for printed model PA-TOOL-01
	IMPLANT LEVEL				
	Conventional impression		Digital impression		
Pop-in transfer Screw included Hexagonal wrench	R PLATFORM 4.8 mm	N PLATFORM 4.0 mm	R PLATFORM 4.8 mm	N PLATFORM 4.0 mm	Digital transfer
	ILT200-R-4 (x4)	ILT200-N-4 (x4)	ILT100-R (x1) ILT100-R-4 (x4)	ILT100-N (x1) ILT100-N-4 (x4)	
Pick-up transfer Non indexed Long and short screws included Hexagonal wrench	ILT100-R (x1) ILT100-R-4 (x4)		ILT100-N (x1) ILT100-N-4 (x4)		156-02-DT-IL
	TA100-R (x1) TA100-R-4 (x4)		TA100-N (x1) TA100-N-4 (x4)		156-02-PA (x1) 156-01-PA (x1)
Abutment analog Hexagonal wrench	TA100-R (x1) TA100-R-4 (x4)		TA100-N (x1) TA100-N-4 (x4)		156-02-PA (x1) 156-01-PA (x1)

MULTI-UNIT ABUTMENTS STERILE 25 N.cm	ABUTMENT	
	Ø 4.8 mm / R	Ø 4.0 mm / N
	TMU0-R	TMU0-N

	MULTI-UNIT				Wrench for analog for printed model PA-TOOL-01
	ABUTMENT LEVEL				
	Conventional impression		Digital impression		
Pop-in transfer Manual	Ø 4.8 mm / R	Ø 4.0 mm / N	Ø 4.8 mm / R	Ø 4.0 mm / N	Digital transfer
	MUT200 (x1) MUT200-4 (x4)	MUNT200 (x1) MUNT200-4 (x4)	151-03-DT-MU	151-04-DT-MUN	
Pick-up transfer Manual Long and short MU screws included	MUT100 (x1) MUT100-4 (x4)		MUNT100 (x1) MUNT100-4 (x4)		151-03-PA (x1) 151-03-PA-4 (x4)
	MUA100 (x1) MUA100-4 (x4)		MUNA100 (x1) MUNA100-4 (x4)		151-04-PA (x1) 151-04-PA-4 (x4)
Analog Manual	MUA100 (x1) MUA100-4 (x4)		MUNA100 (x1) MUNA100-4 (x4)		151-03-PA (x1) 151-03-PA-4 (x4)

	MULTIPLE SCREW RETAINED PROSTHESIS				Wrench for analog for printed model PA-TOOL-01
	IMPLANT LEVEL				
	Conventional impression		Digital impression		
Pop-in transfer Screw included Hexagonal wrench	R PLATFORM 4.8 mm	N PLATFORM 4.0 mm	R PLATFORM 4.8 mm	N PLATFORM 4.0 mm	Digital transfer
	ILT200-R-4 (x4)	ILT200-N-4 (x4)	ILT100-R (x1) ILT100-R-4 (x4)	ILT100-N (x1) ILT100-N-4 (x4)	
Pick-up transfer Non indexed Long and short screws included Hexagonal wrench	ILT100-R (x1) ILT100-R-4 (x4)		ILT100-N (x1) ILT100-N-4 (x4)		156-02-DT-IL
	TA100-R (x1) TA100-R-4 (x4)		TA100-N (x1) TA100-N-4 (x4)		156-02-PA (x1) 156-01-PA (x1)
Abutment analog Hexagonal wrench	TA100-R (x1) TA100-R-4 (x4)		TA100-N (x1) TA100-N-4 (x4)		156-02-PA (x1) 156-01-PA (x1)

MULTI-UNIT Manual	ACCESSORIES & FINAL COPING			
	Ø 4.8 mm / R		Ø 4.0 mm / N	
Protective cap STERILE	MUCAP (x1)	MUCAP-4 (x4)	MUNCAP	
Protective analog	MUA200 (x1)	MUA200-4 (x4)	MUNA200	
Castable coping	MUC300	MUNC300	MUNPAC100	
CoCr castable coping	MUC400	MUNC400	MUNPAC110	
Straight pacific ring			MUNPAC120	
Analog pacific				
Castable pacific				

	TI-BASE			
	SINGLE CEMENTED OR SCREW RETAINED RESTORATIONS			
	FLEXIBASE		CEREC	
Ti-base 25 N.cm MI.6 black screw included Hexagonal wrench	Coronal height: 3.5 mm			
	R PLATFORM 4.8 mm	N PLATFORM 4.0 mm	Ø 5.0 mm Size L Block FOR R PLATFORM 4.8 mm	Ø 4.0 mm Size S Block FOR N PLATFORM 4.0 mm
Cast coping Manual Hexagonal wrench	TFLEX-R	TFLEX-N	TBASEC-R-L	TBASEC-N-S
	TFLEX-R	OPFLEX403		
Laboratory Scan Adapters Manual Hexagonal wrench	156-0X-SAA		Wrench for Scan-Adapter SATOOL-01	

	MULTIPLE SCREW RETAINED RESTORATIONS			
	FLEXIBASE			
	R PLATFORM 4.8 mm		N PLATFORM 4.0 mm	
Ti-base 25 N.cm MI.6 black screw included Hexagonal wrench	Coronal height: 3.5 mm			
	TFLEX-R-P	TFLEX-N-P		
Cast coping Manual Hexagonal wrench	TFLEX-R-P	TFLEX-N-P		
	TFLEX-R-P	TFLEX-N-P		
Laboratory Scan Adapters Manual Hexagonal wrench	156-0X-SAO		Wrench for Scan-Adapter SATOOL-01	

MULTI-UNIT Manual	FLEXIBASE ON MULTI-UNIT			
	Ø 4.8 mm / R		Ø 4.0 mm / N	
Embase 15 N.cm Black screw included MU140	MUFLEX	MUNFLEX		
	MUFLEX	MUNFLEX		
Cast coping Manual Hexagonal wrench	MUFLEX	MUNFLEX		
	MUFLEX	MUNFLEX		
Laboratory Scan Adapters Manual Hexagonal wrench	151-03-SAO		Wrench for Scan-Adapter SATOOL-01	
	151-03-SAO		151-04-SAO	

	TEMPORIZATION			
	SINGLE RESTORATIONS			
	IMPLANT LEVEL			
Temporary abutment indexed 25 N.cm MI.6 titanium screw included Hexagonal wrench	R PLATFORM 4.8 mm		N PLATFORM 4.0 mm	
	TC100-R		TC100-N	

MULTI-UNIT Manual	MULTIPLE SCREW RETAINED PROSTHESIS			
	IMPLANT LEVEL			
Pop-in transfer Screw included Hexagonal wrench	R PLATFORM 4.8 mm	N PLATFORM 4.0 mm	R PLATFORM 4.8 mm	N PLATFORM 4.0 mm
	ILT200-R-4 (x4)	ILT200-N-4 (x4)	ILT100-R (x1) ILT100-R-4 (x4)	ILT100-N (x1) ILT100-N-4 (x4)
Pick-up transfer Non indexed Long and short screws included Hexagonal wrench	ILT100-R (x1) ILT100-R-4 (x4)		ILT100-N (x1) ILT100-N-4 (x4)	
	TA100-R (x1) TA100-R-4 (x4)		TA100-N (x1) TA100-N-4 (x4)	
Abutment analog Hexagonal wrench	TA100-R (x1) TA100-R-4 (x4)		TA100-N (x1) TA100-N-4 (x4)	

MULTI-UNIT Manual	MULTI-UNIT			
	ABUTMENT LEVEL			
Titanium provisional coping 15 N.cm MI.4 MU titanium screw included MU140Z	Ø 4.8 mm / R	Ø 4.0 mm / N	Ø 4.8 mm / R	Ø 4.0 mm / N
	MUC100	MUNC100	MUC200	MUNC200
PEEK provisional coping 15 N.cm MI.4 MU titanium screw included MU140Z	MUC100	MUNC100	MUC200	MUNC200
	MUC200	MUNC200		

INSTRUMENTS & SCREWS		
Black M1.6 prosthetic screw		TS160
Titanium M1.6 prosthetic screw		TS161
Short LAB M1.6 screw		TS163
Long LAB M1.6 screw		TS162
Black M1.6 prosthetic screw (plural)		TS160P
Short LAB M1.6 screw (plural)		TS163P-2 (x2)
Long LAB M1.6 screw (plural)		TS162P-2 (x2)
Short hexagonal wrench		INCHECV
Long hexagonal wrench		INCHELV
XL hexagonal wrench		INCHEXLV
Short hexagonal mandrel		INMHCEV
Long hexagonal mandrel		INMHVELV
XL hexagonal mandrel		INMHXLV
Screwdriver for Scan-Adapters		SATOOL-01
Dynamometrical ratchet wrench		INCCD

Torq Control
Universal dynamometric wrench
15501

Wrench for straight MU abt		MUW100
Mandrel for straight MU abt		MUM100
Long mandrel for straight MU abt		MUM100L
Short prehensile tool for MU		MUWS
M1.6 MU black screw for Axiom® MU abt		OPMU160
M1.6 MU titanium lab screw for Axiom® MU abt (blue)		OPMU161
M1.4 MU black screw for ti co cr prosthesis		MU140
M1.4 MU titanium screw for zirconium prosthesis		MU140Z
M1.4 MU titanium lab screw blue		MU141
MU short lab screw		MUT101
MU long lab screw		MUT102

• MU titanium screw MU140Z included • MU black screw MU140 included • MU laboratory screw MU141 included