

# Dental Wings DWOS Chairside CAD

## v. 2.0.3 - Release Notes



---

---

### IMPROVEMENTS

---

---

Dental Wings is pleased to announce the release of a new software update for DWOS Chairside, the dental CAD/CAM software.

The main improvement in **version 2.0.3** includes an [automatic reconnection](#) to the Dental Wings Intraoral Scanner, [following any network interruption](#). **Download of this newest software version is strongly recommended** if the software is being used in combination with the Dental Wings Intraoral Scanner.

See below for a more detailed summary of latest improvements.

#### Communication with Dental Wings Intraoral Scanner

- Improve the communication between Dental Wings Intraoral Scanner and DWOS Chairside software to ease the reconnection between the 2 software when network interruption occurs.
- Material and implant kit libraries within DWOS Chairside are updated to match the libraries of the Dental Wings Intraoral Scanner.

#### Other

- New classification hierarchy for materials and implant kits belonging to members of the Straumann Group now appear under the same label.

#### BUG FIXES

- Fixed an issue in the *Inbox* which prevented some orders from opening properly.
- The correct information about the blank is displayed when selecting a milling blank.