



Release Notes v. 2023.1¹

We're doing our best to improve Nova every day. See below what we have improved to give you better experience.

New features

[Editable Model ID](#)

You can now change the Model ID to make it easier to keep track of models.

[Visibility options](#)

Models now have visibility options: show/hide and customize transparency for displayed items (model, die, analog or hinges).

[Shortcuts](#)

Show and Hide for upper and lower elements are now available using convenient shortcuts: Alt+C and Alt+V.

[Other language users](#)

With improved translation DWOS is easier to use than ever in German, French, Spanish, Italian, Portuguese.

¹ This version is compatible with Dental Wings standalone software licenses and the following lab scanners: 7Series, 3Series, and Medit T500, T310, T510, and T710.



What we've improved

We've made the following improvements to make your CAD experience better.

New name

Nova now has an official new name: Straumann Nova™.

Faster

Transition time between steps in Model is now faster.

Model only workflow

You can now create a model with analogs and a gingiva mask. You don't need to have a restoration design first.

Functional model

You can create a functional model from a crown or pontic during the design phase, simply by merging the anatomy design with the model.

Windows 11 compatibility

Straumann Nova™ is now compatible with Windows 11.