

Release notes v. 2021.1.1¹

Below you shall find the bug fixes provided in this release.

implant components

Fixed an issue that could potentially led to misaligned components for the following implant platforms

- Straumann RN / WN Implant System
 - Original Straumann RN components
 - Medentika L-Series components
 - Medentika N-Series components
- BEGO Implant System
 - Medentika BS-Series components

Implant libraries

Fixed an issue that prevented access to encrypted implant libraries.

Software update process

Fixed an issue that prevent the update process of Easy while the DWOS server was still running.

Reminder: DWOS & Easy have independent update processes so you must update each application separately once Easy is installed on your system.

¹ This version is compatible with Dental Wings standalone software licenses and the following lab scanners: 7Series, 3Series, and Medit T500, T310, T510, and T710.