

# STRAUMANN UN!Q™

Excel with custom prosthetics.



# SCALE UP WITH STRAUMANN UN!Q™


An outsourced open-platform design and production service for patient-specific prosthetics.

In a time of digital transformation, when resources are scarce and competition is fierce, Straumann UN!Q™ empowers you with a premium service to outsource the planning, design and manufacturing of your custom implant prosthetics on demand, based on your specific needs.

Meeting your growing demand has never been easier with Straumann UN!Q™.



## STRAUMANN UN!Q™ WORKFLOW – OPTIONS FOR YOUR SPECIFIC NEEDS



# EXCEL WITH CUSTOM PROSTHETICS



**Scalable:** grow your business without further investment in CAD resources or CAM equipment.



**Reliable:** deliver consistent, predictable results by leveraging the Straumann® 70-year legacy of Swiss engineering. Over 2 million custom abutments produced in the last 20 years.



**Comprehensive:** extensive portfolio and wide compatibility with over 125 implant platforms and design options to cover every indication.



**Flexible:** pick & choose the on-demand workflow that is tailored to your needs.

Are you a specialist in prosthetic design? Let us support you with ultra-precision manufacturing.

You want to focus your resources on esthetic finishing? Then let us take care of the groundwork to free up your resources.


You have complex edentulous cases to be restored? We can support you with our full-service, premium offering.

# DESIGN & MANUFACTURING SERVICES

Outsourcing your design work is only a few clicks away. Take advantage of our premium on-demand service via Straumann® AXS, the online case ordering and management portal.




Visit Straumann® AXS.



## KEY BENEFITS – STRAUMANN® AXS PORTAL

- Simple and easy to use
- Track & trace your cases all in one place
- Seamless core file workflow
- Quick access to digital training & education opportunities
- Services easily accessible via our e-shop





## THE CORE FILE WORKFLOW

### Order the core file

You can order the design file of your customized abutment or bar.

This service allows you to reduce time by designing and milling the suprastructure (e.g. crown) in-house while waiting for the arrival of your Straumann® order.



After approving the design, you will be able to download the file from the Straumann® AXS Portal.



# CENTRALIZED MANUFACTURING

Design using your preferred CAD system and send directly to Straumann UNIQ™ for production with direct connectivity from 3Shape, exocad and CARES® Visual.

*“Working with your service feels effortless, the flow is very fast and simple. I love the design approval process and the design quality is great.”*

YDL Dental, USA.




Use the Straumann® CAD libraries from the Straumann® Download Center to send directly to our centralized milling facility.



# PORTFOLIO

We offer a comprehensive portfolio for a wide range of indications to restore single, multi and edentulous cases compatible with Straumann® and most third-party implant platforms. Our design service covers the full portfolio to maximize outsourcing opportunities.

Compatible with 125+ implant connections.



## FULL PORTFOLIO RANGE OF REMOVABLE PROSTHETICS

The Straumann UN!Q™ portfolio includes multiple options for removable prosthetic attachment systems including:

- Dolder U Shape & Egg-Shape Bars
- Preci Vertex and Horix Bar
- Ackermann-Bar
- Round Bar
- MP-Clip Bar





## 1 CUSTOM ANGLED SCREW CHANNEL ABUTMENTS

Allows screw channel to be repositioned up to 25° for ideal positioning

- Uncompromised esthetics in anterior region with screw hole placed lingually
- Easy access and improved functionality in the posterior zone

## 2 STRAUMANN® COREBAR™

Straumann® CoreBar™ – milled bar to support zirconia restorations

- Provides larger surface area for bonding vs. titanium-base coping and enhanced distribution of forces throughout the bar
- Lighter than restorations produced with solid zirconia
- Available through design service
- Efficient workflow – download the core file for the bar from the AXS platform and begin milling your zirconia sleeve while the bar is manufactured

## 3 THIMBLE BRIDGE

- Designed for easy maintenance and repair
- Angled screw channels for enhanced esthetics
- Individual crowns allow for a combination of esthetic and high-strength materials

## 4 BASIC & ADVANCED HYBRIDS

- Acrylic hybrids are an ideal solution for ongoing tissue management with easy tissue relining
- Design versatility for retention
- Circumferential ledge provides additional support for teeth to prevent shearing, increasing durability and stability

## 5 NOVALOC® LOCATOR® BAR

- Provide more resistance against mastication forces compared to abutment-supported attachments
- Align attachments into a common insertion axis
- Adjustable height and width
- 0°, 4°, 6° and 8° wall taper

## 6 2 IN 1 TELESCOPIC

- Provides a combination of the passive fit of a fixed restoration and the cleanability of a removable restoration
- Circumferential seal, created between the gingiva and the secondary structure, promotes hygiene and gingival health



# TURNAROUND WORKFLOWS

## DESIGN TURNAROUND TIMES

**SINGLE UNITS  
(ABUTMENTS)  
& 3D MODELS**



**24 hours**

Cut-off time:



**2 pm CET**

**MULTI-UNITS  
(BRIDGES AND BARS)**



**48 hours**

Cut-off time:



**2 pm CET**

## PRODUCTION TURNAROUND TIMES

### CUSTOM ABUTMENT TURNAROUND

(Titanium and CoCr\*)



Case received or design approved prior to 2 pm



Case Shipped

Monday



Tuesday

Tuesday



Wednesday

Wednesday



Thursday

Thursday



Friday

Friday



Monday

### BAR TURNAROUND



Case received or design approved prior to 2 pm



Case Shipped

Monday



Next Wednesday

Tuesday



Next Thursday

Wednesday



Next Friday

Thursday



Next Tuesday

Friday



Next Wednesday

\*Turnaround times for Titanium and CoCr abutments are shown above. An additional 3 days are required for titanium abutments with gold hue.

## RELIABILITY FOR YOUR PEACE OF MIND

- Benefit from state of the art engineered custom prosthetics, manufactured with precision
- Our prosthetics are designed to fit: precision-engineered and milled to perfection. That's why we offer a lifetime warranty for your peace of mind

## AT YOUR SERVICE

- You can rely on our dental technical service for all your questions, they are only one call or click away
- Comprehensive Training & Education service and support can be found conveniently right within the Straumann® AXS order portal



*“The abutments are high-quality and we appreciate the quick turn-around. Placing an order is self-explanatory, we appreciate the flexibility with the workflows and we would definitely recommend the service.”*

**Palma Dental, USA.**

# STRAUMANN UNIQ™ EXCEL WITH CUSTOM PROSTHETICS.



Over 2 million custom abutments  
produced in the last 20 years



70-year legacy of Swiss engineering



Compatible with 125+ implant platforms



Skilled dental technicians at your service,  
only one call/click away



490.898/en/B/00 02/25

## International Headquarters

Institut Straumann AG  
Peter Merian-Weg 12  
CH-4002 Basel, Switzerland  
Phone +41 (0)61 965 11 11  
Fax +41 (0)61 965 11 01  
[www.straumann.com](http://www.straumann.com)

© Institut Straumann AG, 2025. All rights reserved.  
Straumann® and/or other trademarks and logos from Straumann® mentioned herein are  
the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

 **straumann**