

Straumann® Dental Implant System		Ø 3.3 mm		Ø 4.1 mm		Ø 4.8 mm							
System Overview		Standard Plus (SP) NNC Ø 3.5 mm		Standard (S) RN Ø 4.8 mm		Standard (S) RN Ø 4.8 mm		Standard Plus (SP) RN Ø 4.8 mm		Standard (S) WN Ø 6.5 mm		Standard Plus (SP) WN Ø 6.5 mm	
TISSUE LEVEL SURGERY		Standard (S) Mini Implant Ø 3.5 mm		Standard (S) RN Ø 4.8 mm		Standard (S) RN Ø 4.8 mm		Standard Plus (SP) RN Ø 4.8 mm		Standard (S) WN Ø 6.5 mm		Standard Plus (SP) WN Ø 6.5 mm	
Roxolid® SLActive®	4 mm	–	–	–	–	–	033.0435**	–	–	–	–	–	–
	6 mm	–	–	–	–	033.5305	033.5605	–	–	–	–	–	–
	8 mm	–	033.4165	033.5015	033.5115	033.5315	033.5615	–	–	–	–	–	–
	10 mm	–	033.4175	033.5025	033.5125	033.5325	033.5625	–	–	–	–	–	–
	12 mm	–	033.4185	033.5035	033.5135	033.5335	033.5635	–	–	–	–	–	–
	14 mm	–	033.4195	033.5045	033.5145	033.5345	033.5645	–	–	–	–	–	–
	16 mm	–	–	033.5055	–	033.5355	–	–	–	–	–	–	–
	6 mm	–	–	–	–	043.1555	043.1615	–	–	–	–	–	–
	8 mm	–	043.4165	043.2555	043.2605	043.1565	043.1625	–	–	–	–	–	–
	10 mm	042.9445	043.4175	043.2565	043.2615	043.1575	043.1635	–	–	–	–	–	–
	12 mm	042.9455	043.4185	043.2575	043.2625	043.1585	043.1645	–	–	–	–	–	–
	14 mm	042.9465	043.4195	043.2585	043.2635	043.1595	043.1655	–	–	–	–	–	–
	16 mm	–	–	043.2595	–	043.1605	–	–	–	–	–	–	–
Roxolid® SLActive®, guided	8 mm	–	–	033.431G	033.451G	033.531G	033.561G	–	–	–	–	–	–
	10 mm	–	–	033.432G	033.452G	033.532G	033.562G	–	–	–	–	–	–
	12 mm	–	–	033.433G	033.453G	033.533G	033.563G	–	–	–	–	–	–
Round Burs													
Pilot and Twist Drills													
Alignment Pins													
Profile Drills													
Taps													
Closure Caps	NNC		RN		RN				RN		WN		
Healing Caps													



- ABBREVIATIONS/FOOTNOTES**
- 04 Pack of 4
 - AH Abutment height
 - BL Bone Level
 - H Height
 - ND Narrow Diameter
 - NNC Narrow Neck CrossFit®
 - RD Regular Diameter
 - RN Regular Neck
 - S Standard
 - SP Standard Plus
 - V2 Pack of 2
 - V4 Pack of 4
 - WN Wide Neck
- 1) Straumann® CARES® Abutments can be ordered via the Straumann® CARES® Visual Software. Please check portal to see options or at www.cares.straumann.com.
 - 2) Manufactured at the Straumann® CAD/CAM production center
 - 3) Available in 6 different heights
 - 4) Available in 5 different heights
 - 5) Also via Straumann® CARES® Scan & Shape Service
- * Also available without drill stop
- ** For tapping of 4 mm Standard Plus Short Ø 4.1 mm implant use BL/TE Tap 026.4310
- *** For tapping of 4 mm Standard Plus Short Ø 4.8 mm implant use BL/TE Tap 026.6310
- Crown
 - Bridge
 - Bar
- Novaloc® is a registered trademark of Valoc AG, Switzerland.
- Product images are not to scale.

International Headquarters
 Institut Straumann AG
 Peter Merian-Weg 12
 CH-4002 Basel, Switzerland
 Phone +41 (0)61 965 11 11
 Fax +41 (0)61 965 11 01
www.straumann.com

© Institut Straumann AG, 2023. All rights reserved.
 Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates. All rights reserved.

