

SELECTION GUIDE

Straumann® Modular Cassette




CONTENT

A Module Full List	3
Standard Plus (SP), Standard Plus Short (SPS), Bone Level (BL)	8
Standard Plus (SP) guided, Standard Plus Short (SPS) guided, Bone Level (BL) guided	12
Straumann® PURE Ceramic Implant System (CIM/CI)	20
Bone Level Tapered (BLT)	26
Bone Level Tapered (BLT) guided	30
Tissue Level X (TLX), Bone level X (BLX)	38
Tissue Level X (TLX) guided, Bone Level X (BLX) guided	42

PUPOSE OF THIS GUIDE

The selection guide helps to define the required Modular Cassette parts and instruments needed for a specific indication. Each chapter includes an order list as well as a set-up overview which shows how to assemble the parts.




* For Straumann® PURE Ceramic Implant System using guided surgery, please use the SP/SPS/BL guided surgery instrument setup (page 12) for guided implant bed preparation. For implant placement use the CIM/CI freehand instrument set-up (page 20).

A MODULE FULL LIST


This chapter shows all mandatory *and* optional instruments for the A Module.






A MODULE




A Module		
Art. No.	Image	Product
041.761		Straumann® Modular Cassette, A Module

RATCHET TRAY




Ratchet Tray		
Art. No.	Image	Product
041.766		A Modul Ratchet Tray
046.119		Ratchet
046.049		Torque Control Device, 15 Ncm, 35 Ncm
066.1100		Torque Control Device, 35 Ncm, 50 Ncm, 80 Ncm
046.064		Holding Key

RELEASE AID FOR LOXIM® TRAY




Release Aid Tray			
Art. No.	Image	Product	Indication
041.767		Release Aid Tray	
026.0073		Release Aid S for Loxim®	BLT SC
026.2558		Release Aid N for Loxim®	S, SP, TE, BL, BLT
026.4558		Release Aid R/W for Loxim®	S, SP, TE, BL, BLT

STOP KEY TRAY




Stop Key Tray		
Art. No.	Image	Product
041.768		Stop Key Tray
034.006		Stop Key, guided





TWEEZER TRAY




Tweezer Tray		
Art. No.	Image	Product
041.770		Tweezer Tray
046.095		Tweezers with titanium plasma-coated tips



PLANNING AID TRAY



Planning Aid Tray			
Art. No.	Image	Product	Indication
041.769		Planning Aid Tray	
046.147		Diagnostic T	S, SP, TE, CIM/CI, BL, BLT
046.148		Tissue Level Implant Distance Indicator	S, SP, TE, CIM/CI
026.0901		Bone Level Implant Distance Indicator	BL, BLT

IMPLANT DEPTH GAUGE TRAY



Implant Depth Gauge Tray		
Art. No.	Image	Product
041.771		Implant Depth Gauge Tray
066.2000		Implant Depth Gauge

SCS SCREWDRIVER SET



SCS Screwdriver Set		
Art. No.	Image	Product
041.764		Grommet Tray, 3 small + 3 large
046.400		SCS Screwdriver for Ratchet, extra-short
046.401		SCS Screwdriver for Ratchet, short
046.402		SCS Screwdriver for Ratchet, long
046.410		SCS Screwdriver for Handpiece, extra-short
046.411		SCS Screwdriver for Handpiece, short
046.412		SCS Screwdriver for Handpiece, long

AS SCREWDRIVER SET



AS Screwdriver Set		
Art. No.	Image	Product
041.764		Grommet Tray, 3 small + 3 large
046.789		AS Screwdriver for Handpiece, extra short
046.790		AS Screwdriver for Handpiece, short
046.791		AS Screwdriver for Handpiece, long
046.786		AS Screwdriver for Ratchet, extra-short*
046.787		AS Screwdriver for Ratchet, short*
046.788		AS Screwdriver for Ratchet, long*

* Store the handling aid of the AS Screwdriver for Ratchet separately in the 041.765 Grommet tray

BONE PROFILER SET



Bone Profiler Set			
Art. No.	Image	Product	Indication
041.762		Grommet Tray 6 small	
044.304		Bone Profiler RN/RT	S, SP, TLX
044.305		Bone Profiler WN/WT	
026.0070		Bone Profiler, oval Platform SC	BLT SC Ø2.9 mm
026.0022		BL Bone Profiler 1	BL, BLT, BLX
026.0023		BL Bone Profiler 2	
026.0024		BL Bone Profiler 3	

48H TOOL



48h Tool			
Art. No.	Image	Product	Indication
041.763		Grommet Tray 6 large	
80103		48 hour Explantation Device RN/WN	S, SP, TE
80087		48 hour Explantation Device RN/WN	
031.080		48 hour Explantation Device, Straumann® PURE Ceramic Implant (Monotype) RD	CIM/CI
031.081		48 hour Explantation Device, Straumann® PURE Ceramic Implant (Monotype) ND	
032.053		48h Explantation Device, Straumann® PURE Ceramic Implant RD	
026.0072		48 hour Explantation Device SC	BLT
026.2048		48 hour Explantation Device NC	BL, BLT
026.4048		48 hour Explantation Device RC	
049.361		48 hour Explantation Device NNC	NNC
066.4205		TorcFit™ Implant Driver for Ratchet, short	BLX, TLX
066.4206		TorcFit™ Implant Driver for Ratchet, long	

* Note: CIM/CI PURE 48h explantation devices are for single use only.



Order List: Standard Plus, Standard Plus Short, Bone Level

A MODULE – ORDER LIST

SP Standard Plus / SPS Standard Plus Short / BL Bone Level		
Art. No.	Image	Product
041.761		Straumann® Modular Cassette, A Module
041.766		A Modul Ratchet Tray
046.119		Ratchet
046.049		Torque Control Device
046.064		Holding Key
041.767		Release Aid Tray
026.2558		Release Aid N for Loxim®
026.4558		Release Aid R/W for Loxim®
041.764		Grommet Tray, 3 small + 3 large
046.400		SCS Screwdriver for Ratchet, extra-short
046.401		SCS Screwdriver for Ratchet, short
046.402		SCS Screwdriver for Ratchet, long
046.410		SCS Screwdriver for Handpiece, extra-short
046.411		SCS Screwdriver for Handpiece, short
046.412		SCS Screwdriver for Handpiece, long

A MODULE TISSUE LEVEL SPECIFIC TOOLS

Add for Standard Plus and Standard Plus Short Implant line

SP Standard Plus / SPS Standard Plus Short / BL Bone Level		
Art. No.	Image	Product
041.764		Grommet Tray, 3 small + 3 large
044.304		Bone Profiler RN/RT
044.305		Bone Profiler WN/WT
80103		48 hour Explantation Device RN/WN
80087		48 hour Explantation Device RN/WN

A MODULE BONE LEVEL SPECIFIC TOOLS

Add for Bone Level Implant Line

SP Standard Plus / SPS Standard Plus Short / BL Bone Level		
Art. No.	Image	Product
041.764		Grommet Tray, 3 small + 3 large
026.0022		BL Bone Profiler 1
026.0023		BL Bone Profiler 2
026.0024		BL Bone Profiler 3
026.2048		48 hour Explantation Device NC
026.4048		48 hour Explantation Device RC

For additional tools (planning aid, tweezers, AS Screwdriver), please refer to page 3 *A Module full list*.



Order List: Standard Plus, Standard Plus Short, Bone Level

B MODULE - ORDER LIST

SP Standard Plus / SPS Standard Plus Short / BL Bone Level		
Art. No.	Image	Product
041.776		Straumann® Modular Cassette, B Module
041.779		B Module, Parallel Walled Tray
040.563		Drill extender
044.003		Round bur, Ø 2.3 mm
044.004		Round bur, Ø 3.1 mm
026.0056		Needle drill, long, Ø 1.6 mm
044.210		Pilot Drill 1, short, Ø 2.2 mm
044.214		Pilot Drill 2, short, Ø 2.8 mm
044.250		Twist Drill PRO, short, Ø 3.5 mm
044.254		Twist Drill PRO, short, Ø 4.2 mm
044.211		Pilot Drill 1, long, Ø 2.2 mm
044.215		Pilot Drill 2, long, Ø 2.8 mm
044.251		Twist Drill PRO, long, Ø 3.5 mm
044.255		Twist Drill PRO, long, Ø 4.2 mm
046.799		Alignment Pin, Ø 2.2 mm
046.800		Depth Gauge, Ø 2.8 mm
046.802		Depth Gauge, Ø 3.5 mm
046.804		Depth Gauge, Ø 4.2 mm

046.460		Adapter for Ratchet, extra-short
046.461		Adapter for Ratchet, short
046.470		Adapter for Handpiece, extra-short
046.471		Adapter for Handpiece, short

B MODULE TISSUE LEVEL SPECIFIC TOOLS

Add for Standard Plus and Standard Plus Short Implant line

SP Standard Plus / SPS Standard Plus Short / BL Bone Level		
Art. No.	Image	Product
034.332		SP Profile Drill, short, Ø 3.3 mm, FIBA compatible
034.333		SP Profile Drill, short, Ø 4.1 mm, FIBA compatible
044.084		SP Profile Drill, short, Ø 4.8 mm, FIBA compatible
034.351		S/SP Tap, short, Ø 3.3 mm, FIBA compatible
034.352		S/SP Tap, short, Ø 4.1 mm, FIBA compatible
034.353		S/SP Tap, short, Ø 4.8 mm, FIBA compatible

B MODULE BONE LEVEL SPECIFIC TOOLS


Add for Bone Level Implant Line

SP Standard Plus / SPS Standard Plus Short / BL Bone Level		
Art. No.	Image	Product
034.327		BL/NNC Profile Drill, short, Ø 3.3 mm, FIBA compatible
034.328		BL Profile Drill, short, Ø 4.1 mm, FIBA compatible
034.329		BL Profile Drill, short, Ø 4.8 mm, FIBA compatible
034.348		BL/NNC Tap, short, Ø 3.3 mm, FIBA compatible
034.349		BL Tap, short, Ø 4.1 mm, FIBA compatible
034.350		BL Tap, short, Ø 4.8 mm, FIBA compatible



A Module Set-up: Standard Plus, Standard Plus Short, Bone Level


A Module
041.761






B Module Set-up: Standard Plus, Standard Plus Short, Bone Level






Parallel Walled Tray
041.779



Bone Level Tools



Tissue Level Tools





Order List: Standard Plus Guided, Standard Plus Short Guided, Bone Level Guided

A MODULE – ORDER LIST

SP Standard Plus / SPS Standard Plus Short / BL Bone Level		
Art. No.	Image	Product
041.761		Straumann® Modular Cassette, A Module
041.766		A Modul Ratchet Tray
046.119		Ratchet
046.049		Torque Control Device
046.064		Holding Key
041.767		Release Aid Tray
026.2558		Release Aid N for Loxim®
026.4558		Release Aid R/W for Loxim®
041.764		Grommet Tray, 3 small + 3 large
046.400		SCS Screwdriver for Ratchet, extra-short
046.401		SCS Screwdriver for Ratchet, short
046.402		SCS Screwdriver for Ratchet, long
046.410		SCS Screwdriver for Handpiece, extra-short
046.411		SCS Screwdriver for Handpiece, short
046.412		SCS Screwdriver for Handpiece, long

A MODULE TISSUE LEVEL SPECIFIC TOOLS

Add for Standard Plus and Standard Plus Short Implant line

SP Standard Plus / SPS Standard Plus Short / BL Bone Level		
Art. No.	Image	Product
041.764		Grommet Tray, 3 small + 3 large
044.304		Bone Profiler RN/RT
044.305		Bone Profiler WN/WT
80103		48 hour Explantation Device RN/WN
80087		48 hour Explantation Device RN/WN

A MODULE BONE LEVEL SPECIFIC TOOLS

Add for Bone Level Implant Line

SP Standard Plus / SPS Standard Plus Short / BL Bone Level		
Art. No.	Image	Product
041.764		Grommet Tray, 3 small + 3 large
026.0022		BL Bone Profiler 1
026.0023		BL Bone Profiler 2
026.0024		BL Bone Profiler 3
026.2048		48 hour Explantation Device NC
026.4048		48 hour Explantation Device RC

For additional tools (planning aid, tweezers, AS Screwdriver), please refer to page 3 *A Module full list*.



Order List: Standard Plus Guided, Standard Plus Short Guided, Bone Level Guided

B MODULE - ORDER LIST

SP Standard Plus / SPS Standard Plus Short / BL Bone Level		
Art. No.	Image	Product
041.776		Straumann® Modular Cassette, B Module
041.782		B Module, Guided Parallel Walled Tray
034.010		Mucosa Punch, Ø 3.4 mm
034.011		Mucosa Punch, Ø 4.0 mm
034.012		Mucosa Punch, Ø 4.7 mm
034.215		Milling Cutter, Ø 2.8 mm
034.415		Milling Cutter, Ø 3.5 mm
034.615		Milling Cutter, Ø 4.2 mm
034.123		Guided Pilot Drill, short, Ø 2.2 mm
034.223		Guided Twist drill PRO, short, Ø 2.8 mm
034.423		Guided Twist drill PRO, short, Ø 3.5 mm
034.623		Guided Twist drill PRO, short, Ø 4.2 mm
034.126		Guided Pilot Drill, medium, Ø 2.2 mm
034.226		Guided Twist drill PRO, medium, Ø 2.8 mm
034.426		Guided Twist drill PRO, medium, Ø 3.5 mm
034.626		Guided Twist drill PRO, medium, Ø 4.2 mm
034.129		Guided Pilot Drill, long, Ø 2.2 mm
034.229		Guided Twist drill PRO, long, Ø 2.8 mm
034.429		Guided Twist drill PRO, long, Ø 3.5 mm
034.629		Guided Twist drill PRO, long, Ø 4.2 mm
046.799		Alignment Pin, Ø 2.2 mm, BLT
046.800		Depth Gauge, Ø 2.8 mm, BLT
046.802		Depth Gauge, Ø 3.5 mm, BLT
046.804		Depth Gauge, Ø 4.2 mm, BLT
034.354		Guided Adapter for profile drill, FIBA, H2
034.355		Guided Adapter for profile drill, FIBA, H4
034.356		Guided Adapter for profile drill, FIBA, H6

034.357		Guided Adapter for tap, FIBA, H2
034.358		Guided Adapter for tap, FIBA, H4
034.359		Guided Adapter for tap, FIBA, H6

B MODULE BONE LEVEL SPECIFIC TOOLS

Add for Bone Level Implant Line

SP Standard Plus / SPS Standard Plus Short / BL Bone Level		
Art. No.	Image	Product
034.327		BL/NNC Profile Drill, short, Ø 3.3 mm, FIBA compatible
034.328		BL Profile Drill, short, Ø 4.1 mm, FIBA compatible
034.329		BL Profile Drill, short, Ø 4.8 mm, FIBA compatible
034.348		BL/NNC Tap, short, Ø 3.3 mm, FIBA compatible
034.349		BL Tap, short, Ø 4.1 mm, FIBA compatible
034.350		BL Tap, short, Ø 4.8 mm, FIBA compatible
026.0083		BL/BLT Guided Adapter, for handpiece
026.0084		BL/BLT Guided Adapter, for ratchet

B MODULE TISSUE LEVEL SPECIFIC TOOLS












Add for Standard Plus and Standard Plus Short Implant line

SP Standard Plus / SPS Standard Plus Short / BL Bone Level		
Art. No.	Image	Product
034.332		SP Profile Drill, short, Ø 3.3 mm, FIBA compatible
034.333		SP Profile Drill, short, Ø 4.1 mm, FIBA compatible
044.084		SP Profile Drill, short, Ø 4.8 mm, FIBA compatible
034.351		S/SP Tap, short, Ø 3.3 mm, FIBA compatible
034.352		S/SP Tap, short, Ø 4.1 mm, FIBA compatible
034.353		S/SP Tap, short, Ø 4.8 mm, FIBA compatible
046.708		SP/NNC Guided Adapter, for handpiece
046.710		SP/NNC Guided Adapter, for ratchet



Order List: Standard Plus Guided, Standard Plus Short Guided, Bone Level Guided

C MODULE – ORDER LIST


TLX/BLX Guided		
Art. No.	Image	Product
041.772		Straumann® Modular Cassette, C Module, Guided surgery
026.0147		Drill handle, Ø 2.2 mm, 1 mm/3 mm
026.0148		Drill handle, Ø 2.8 mm, 1 mm/3 mm
026.0150		Drill handle, Ø 3.5 mm, 1 mm/3 mm
026.0152		Drill handle, Ø 4.2 mm, 1 mm/3 mm
034.284		Drill for Template Fixation Pin, Ø 1.3 mm
034.282		Template Fixation Pin, Ø 1.3 mm
034.298		Template Fixation Pin, Ø 2.8/2.8 mm
034.285		Template Fixation Pin, Ø 5/2.8 mm
034.287		Template Fixation Pin, Ø 5/3.5 mm
034.289		Template Fixation Pin, Ø 5/4.2 mm





A Module Set-up: Standard Plus ^{1/3} Guided, Standard Plus Short Guided, Bone Level Guided

A Module
 041.761






B Module Set-up: Standard Plus ^{2/3} Guided, Standard Plus Short Guided, Bone Level Guided

Guided Parallel Walled Tray 041.782




Guided Adapter for profile drill
034.354
034.355
034.356



Guided Adapter for tap
034.357
034.358
034.359



Bone Level Tools




Tissue Level Tools





A Module Set-up: Standard Plus ^{3/3} Guided, Standard Plus Short Guided, Bone Level Guided

C Module Guided Surgery
 041.772







Order List: Straumann® PURE Ceramic Implant System

A MODULE – ORDER LIST

Straumann® PURE Ceramic Implant System		
Art. No.	Image	Product
041.761		Straumann® Modular Cassette, A Module
041.766		A Modul Ratchet Tray
046.119		Ratchet
046.049		Torque Control Device
046.064		Holding Key
041.764		Grommet Tray, 3 small + 3 large
046.400		SCS Screwdriver for Ratchet, extra-short
046.401		SCS Screwdriver for Ratchet, short
046.402		SCS Screwdriver for Ratchet, long
046.410		SCS Screwdriver for Handpiece, extra short
046.411		SCS Screwdriver for Handpiece, short
046.412		SCS Screwdriver for Handpiece, long
041.763		Grommet Tray, 6 large
031.080		48 hour Explantation Device, Straumann® PURE Ceramic Implant (Monotype) RD*
031.081		48 hour Explantation Device, Straumann® PURE Ceramic Implant (Monotype) ND*
032.053		48h Explantation Device, Straumann® PURE Ceramic Implant RD*

* Note: PURE 48h explantation device are for single use only.

For additional tools (planning aid, tweezers, AS Screwdriver), please refer to page 3 *A Module full list*.



Order List: Straumann® PURE Ceramic Implant System

B MODULE - ORDER LIST

Straumann® PURE Ceramic Implant System		
Art. No.	Image	Product
041.776		Straumann® Modular Cassette, B Module, Base + Lid
041.778		B Module, Basic Tray
040.563		Drill extender
044.003		Round bur, Ø 2.3 mm
044.004		Round bur, Ø 3.1 mm
026.0056		Needle drill long, Ø 1.6 mm
044.210		Pilot Drill 1, short, Ø 2.2 mm
044.214		Pilot Drill 2, short, Ø 2.8 mm
044.250		Twist Drill PRO, short, Ø 3.5 mm
044.254		Twist Drill PRO, short, Ø 4.2 mm
044.211		Pilot Drill 1, long, Ø 2.2 mm
044.215		Pilot Drill 2, long, Ø 2.8 mm
044.251		Twist Drill PRO, long, Ø 3.5 mm
044.255		Twist Drill PRO, long, Ø 4.2 mm
034.327		BL/NNC Profile Drill, short, Ø 3.3 mm, FIBA compatible
034.328		BL Profile Drill, short, Ø 4.1 mm, FIBA compatible
034.329		BL Profile Drill, short, Ø 4.8 mm, FIBA compatible
034.348		BL/NNC Tap, short, Ø 3.3 mm, FIBA compatible
034.349		BL Tap, short, Ø 4.1 mm, FIBA compatible
034.350		BL Tap, short, Ø 4.8 mm, FIBA compatible

046.799		Alignment Pin, Ø 2.2 mm, BLT
046.800		Depth Gauge Ø 2.8 mm, BLT
046.802		Depth Gauge, Ø 3.5 mm, BLT
046.804		Depth Gauge, Ø 4.2 mm, BLT
046.460		Adapter for Ratchet, extra-short
046.461		Adapter for Ratchet, short
046.470		Adapter for Handpiece, extra-short
046.471		Adapter for Handpiece, short


C MODULE - ORDER LIST

Straumann® PURE Ceramic Implant System		
Art. No.	Image	Product
041.773		Straumann® Modular, Cassette, C Module, Spare Tools
031.102		Position Indicator, PURE ND, AH 4 mm, Ø 2.2 mm
031.103		Position Indicator, PURE ND, AH 4 mm, Ø 2.8 mm
031.112		Position Indicator, PURE ND, AH 5.5 mm, Ø 2.2 mm
031.113		Position Indicator, PURE ND, AH 5.5 mm, Ø 2.8 mm
031.123		Position Indicator, PURE RD, AH 4 mm, Ø 2.2 mm
031.125		Position Indicator, PURE RD, AH 4 mm, Ø 3.5 mm
031.143		Position Indicator, PURE RD, AH 5.5 mm, Ø 2.2 mm
031.145		Position Indicator, PURE RD, AH 5.5 mm, Ø 3.5 mm



A Module Set-up: Straumann® PURE Ceramic Implant System


A Module
041.761





B Module Set-up: Straumann® PURE Ceramic Implant System


Basic Tray 041.778





C Module Set-up: Straumann® PURE Ceramic Implant System

C Module Spare Tools
041.773



Position Indicator	031.103	031.125
031.102	031.113	031.145
031.112		
031.123		
031.143		



BLT Order List: Bone Level Tapered

A MODULE – ORDER LIST

BLT Bone Level Tapered		
Art. No.	Image	Product
041.761		Straumann® Modular Cassette, A Module
041.763		Grommet Tray, 6 large
026.2048		48 hour tool NC, BLT
026.4048		48 hour tool RC, BLT
041.766		A Modul Ratchet Tray
046.119		Ratchet
046.049		Torque Control Device
046.064		Holding Key
041.767		Release Aid Tray
026.2558		Release Aid N for Loxim®
026.4558		Release Aid R/W for Loxim®
041.764		Grommet Tray, 3 small + 3 large
046.400		SCS Screwdriver for Ratchet, extra-short
046.401		SCS Screwdriver for Ratchet, short
046.402		SCS Screwdriver for Ratchet, long
046.410		SCS Screwdriver for Handpiece, extra-short
046.411		SCS Screwdriver for Handpiece, short
046.412		SCS Screwdriver for Handpiece, long

041.762		Grommet Tray, 6 small
026.0022		BL Bone Profiler 1
026.0023		BL Bone Profiler 2
026.0024		BL Bone Profiler 3

A MODULE BLT Ø 2.9 SPECIFIC TOOLS

Add for BLT Ø 2.9 implant line








BLT Bone Level Tapered Ø 2.9		
Art. No.	Image	Product
041.767		Release Aid Tray*
026.0073		Release Aid S for Loxim®
041.764		Grommet Tray, 3 small + 3 large*
026.0070		Bone Profiler, oval Platform, BLT SC*
026.0072		48 hour tool SC, BLT SC*
025.0042		SC Adapter for Handpiece, long, BLT SC*
025.0043		SC Adapter for Ratchet, long, BLT SC*













*For Bone Level Tapered Small CrossFit® implant line only

For additional tools (planning aid, tweezers, AS Screwdriver), please refer to page 3 *A Module full list*.

BLT Order List: Bone Level Tapered

B MODULE - ORDER LIST

BLT Bone Level Tapered		
Art. No.	Image	Product
041.776		Straumann® Modular Cassette, B Module, Base + Lid
041.778		B Module, Basic Tray
040.563		Drill extender
044.003		Round bur, Ø2.3 mm
044.004		Round bur, Ø3.1 mm
026.0056		Needle drill, long, Ø1.6 mm
066.1301		X VeloDrill™, short, Ø2.2 mm
066.1302		X VeloDrill™, short, Ø2.8 mm
066.1304		X VeloDrill™, short, Ø3.5 mm
066.1306		X VeloDrill™, short, Ø4.2 mm
066.1701		X VeloDrill™, long, Ø2.2 mm
066.1702		X VeloDrill™, long, Ø2.8 mm
066.1704		X VeloDrill™, long, Ø3.5 mm
066.1706		X VeloDrill™, long, Ø4.2 mm
026.0061		BLT Profile Drill, short, Ø2.9 mm
034.324		BLT Profile Drill, short, Ø3.3 mm, FIBA compatible
034.325		BLT Profile Drill, short, Ø4.1 mm, FIBA compatible
034.326		BLT Profile Drill, short, Ø4.8 mm, FIBA compatible

026.0063		BLT Tap, long, Ø2.9 mm
034.345		BLT Tap, short, Ø3.3 mm, FIBA compatible
034.346		BLT Tap, short, Ø4.1 mm, FIBA compatible
034.347		BLT Tap, short, Ø4.8 mm, FIBA compatible
046.799		Alignment Pin, Ø2.2 mm
046.800		Depth Gauge, Ø2.8 mm
046.802		Depth Gauge, Ø3.5 mm
046.804		Depth Gauge, Ø4.2 mm
046.460		Adapter for Ratchet, extra-short
046.461		Adapter for Ratchet, short
046.470		Adapter for Handpiece, extra-short
046.471		Adapter for Handpiece, short

BLT A Module Set-up: Bone Level Tapered


A Module
041.761

Ratchet Tray
041.766

- Ratchet 046.119
- Torque Control 046.049
- Holding Key 046.064

Release Aid Tray
041.767

- Release Aid N 026.2558
- Release Aid R/W 026.4558
- Release Aid S* 026.0073



SC Adapter*
for Ratchet
025.0043

48 hour Tool SC*
026.0072

- SCS Screwdriver for ratchet 046.400, 046.401, 046.402
- SCS Screwdriver for Handpiece 046.410, 046.411, 046.412

SCS Screwdriver Set
Grommet Tray 3 small + 3 large
041.764

- BL Bone Profiler 026.0022, 026.0023, 026.0024

Bone Profiler Set
Grommet Tray 6 small
041.762

- 48 hour Tool NC/RC 026.2048, 026.4048


48h tool Set
Grommet Tray 6 large
041.763

- SC Adapter* for Handpiece 025.0042
- Bone Profiler SC* 026.0070

SC Auxiliary Set*
Grommet Tray 3 small + 3 large
041.764

BLT B Module Set-up: Bone Level Tapered





Basic Tray
041.778



BLT Order List: Bone Level Tapered Guided

A MODULE – ORDER LIST

BLT Bone Level Tapered		
Art. No.	Image	Product
041.761		Straumann® Modular Cassette, A Module
041.763		Grommet Tray, 6 large
026.2048		48 hour tool NC, BLT
026.4048		48 hour tool RC, BLT
041.766		A Modul Ratchet Tray
046.119		Ratchet
046.049		Torque Control Device
046.064		Holding Key
041.767		Release Aid Tray
026.2558		Release Aid N for Loxim®
026.4558		Release Aid R/W for Loxim®
041.764		Grommet Tray, 3 small + 3 large
046.400		SCS Screwdriver for Ratchet, extra-short
046.401		SCS Screwdriver for Ratchet, short
046.402		SCS Screwdriver for Ratchet, long
046.410		SCS Screwdriver for Handpiece, extra-short
046.411		SCS Screwdriver for Handpiece, short
046.412		SCS Screwdriver for Handpiece, long

041.762		Grommet Tray, 6 small
026.0022		BL Bone Profiler 1
026.0023		BL Bone Profiler 2
026.0024		BL Bone Profiler 3

For additional tools (planning aid, tweezers, AS Screwdriver), please refer to page 3 *A Module full list*.












BLT Order List: Bone Level Tapered Guided

B MODULE - ORDER LIST

BLT Bone Level Tapered		
Art. No.	Image	Product
041.776		Straumann® Modular Cassette, B Module, Base + Lid
041.781		B Module, Guided Basic Tray
034.010		Mucosa Punch, Ø 3.4 mm
034.011		Mucosa Punch, Ø 4.0 mm
034.012		Mucosa Punch, Ø 4.7 mm
034.215		Milling Cutter, Ø 2.8 mm
034.415		Milling Cutter, Ø 3.5 mm
034.615		Milling Cutter, Ø 4.2 mm
066.1301		X VeloDrill™, short, Ø 2.2 mm
066.1302		X VeloDrill™, short, Ø 2.8 mm
066.1304		X VeloDrill™, short, Ø 3.5 mm
066.1306		X VeloDrill™, short, Ø 4.2 mm
066.1501		X VeloDrill™, medium, Ø 2.2 mm
066.1502		X VeloDrill™, medium, Ø 2.8 mm
066.1504		X VeloDrill™, medium, Ø 3.5 mm
066.1506		X VeloDrill™, medium, Ø 4.2 mm
066.1701		X VeloDrill™, long, Ø 2.2 mm
066.1702		X VeloDrill™, long, Ø 2.8 mm
066.1704		X VeloDrill™, long, Ø 3.5 mm
066.1706		X VeloDrill™, long, Ø 4.2 mm

034.354		Guided Adapter for profile drill, FIBA, H2
034.355		Guided Adapter for profile drill, FIBA, H4
034.356		Guided Adapter for profile drill, FIBA, H6
034.357		Guided Adapter for tap, H2
034.358		Guided Adapter for tap, H4
034.359		Guided Adapter for tap, H6
034.324		BLT Profile Drill, short, Ø 3.3 mm, FIBA compatible
034.325		BLT Profile Drill, short, Ø 4.1 mm, FIBA compatible
034.326		BLT Profile Drill, short, Ø 4.8 mm, FIBA compatible
034.345		BLT Tap, short, Ø 3.3 mm, FIBA compatible
034.346		BLT Tap, short, Ø 4.1 mm, FIBA compatible
034.347		BLT Tap, short, Ø 4.8 mm, FIBA compatible
046.799		Alignment Pin, Ø 2.2 mm
046.800		Depth Gauge, Ø 2.8 mm
046.802		Depth Gauge, Ø 3.5 mm
046.804		Depth Gauge, Ø 4.2 mm
026.0083		BL/BLT Guided Adapter, for handpiece
026.0084		BL/BLT Guided Adapter, for ratchet


C MODULE – ORDER LIST

TLX/BLX Guided		
Art. No.	Image	Product
041.772		Straumann® Modular Cassette, C Module, Guided surgery
026.0147		Drill handle, Ø 2.2 mm, 1 mm/3 mm
026.0148		Drill handle, Ø 2.8 mm, 1 mm/3 mm
026.0150		Drill handle, Ø 3.5 mm, 1 mm/3 mm
026.0152		Drill handle, Ø 4.2 mm, 1 mm/3 mm
034.284		Drill for Template Fixation Pin, Ø 1.3 mm
034.282		Template Fixation Pin, Ø 1.3 mm
034.298		Template Fixation Pin, Ø 2.8/2.8 mm
034.285		Template Fixation Pin, Ø 5/2.8 mm
034.287		Template Fixation Pin, Ø 5/3.5 mm
034.289		Template Fixation Pin, Ø 5/4.2 mm




BLT A Module Set-up: Bone Level Tapered Guided

A Module
041.761




BLT B Module Set-up: Bone Level Tapered Guided

Guided Basic Tray 041.781



BLT C Module Set-up: Bone Level Tapered Guided

C Module Guided Surgery
041.772





A MODULE – ORDER LIST

TLX/BLX		
Art. No.	Image	Product
041.761		Straumann® Modular Cassette, A Module
041.766		A Modul Ratchet Tray
046.119		Ratchet
066.1100		Torque Control Device
046.064		Holding Key
041.764		Grommet Tray, 3 small + 3 large
046.400		SCS Screwdriver for Ratchet, extra-short
046.401		SCS Screwdriver for Ratchet, short
046.402		SCS Screwdriver for Ratchet, long
046.410		SCS Screwdriver for Handpiece, extra-short
046.411		SCS Screwdriver for Handpiece, short
046.412		SCS Screwdriver for Handpiece, long
041.771		Implant Depth Gauge Tray
066.2000		Implant Depth Gauge

041.762		Grommet Tray, 6 small
026.0022		BL Bone Profiler 1
026.0023		BL Bone Profiler 2
026.0024		BL Bone Profiler 3
044.304		Bone Profiler RN/RT
044.305		Bone Profiler WN/WT
041.762		A Module, Grommet Tray, 6 small
036.3300		RT Profile Drill, short, for implants Ø 3.75/Ø 4.8 mm
036.3301		RT Profile Drill, long, for implants Ø 3.75/Ø 4.8 mm
036.3302		WT Profile Drill, short, for implants Ø 5.5/Ø 6.5 mm
036.3303		WT Profile Drill, long, for implants Ø 5.5/Ø 6.5 mm

For additional tools, tweezers, AS Screwdriver, please refer to page 3 *A Module full list*.

TLX BLX Order List: Tissue Level X, Bone Level X


B MODULE - ORDER LIST

TLX/BLX		
Art. No.	Image	Product
041.776		Straumann® Modular Cassette, B Module, Base + Lid
041.777		B Module, Fully Tapered Tray
044.003		Round bur, Ø 2.3 mm
044.004		Round bur, Ø 3.1 mm
026.0056		Needle drill, long, Ø 1.6 mm
066.1301		X VeloDrill™, short, Ø 2.2 mm
066.1302		X VeloDrill™, short, Ø 2.8 mm
066.1303		X VeloDrill™, short, Ø 3.2 mm
066.1304		X VeloDrill™, short, Ø 3.5 mm
066.1305		X VeloDrill™, short, Ø 3.7 mm
066.1306		X VeloDrill™, short, Ø 4.2 mm
066.1307		X VeloDrill™, short, Ø 4.7 mm
066.1308		X VeloDrill™, short, Ø 5.2 mm
066.1309		X VeloDrill™, short, Ø 6.2 mm
066.1701		X VeloDrill™, long, Ø 2.2 mm
066.1702		X VeloDrill™, long, Ø 2.8 mm
066.1703		X VeloDrill™, long, Ø 3.2 mm
066.1704		X VeloDrill™, long, Ø 3.5 mm
066.1705		X VeloDrill™, long, Ø 3.7 mm
066.1706		X VeloDrill™, long, Ø 4.2 mm, BLX
066.1707		X VeloDrill™, long, Ø 4.7 mm, BLX

046.799		Alignment Pin, Ø 2.2 mm
046.800		Depth Gauge, Ø 2.8 mm
046.801		Depth Gauge, Ø 3.2 mm
046.802		Depth Gauge, Ø 3.5 mm
046.803		Depth Gauge, Ø 3.7 mm
046.804		Depth Gauge, Ø 4.2 mm
046.805		Depth Gauge, Ø 4.7 mm
046.806		Depth Gauge, Ø 5.2 mm
046.807		Depth Gauge, Ø 6.2 mm
066.4201		Implant driver for Ratchet, short
066.4207		Implant driver for Ratchet, medium
066.4202		Implant driver for Ratchet, long
066.4101		Implant driver for Handpiece, short
066.4107		Implant driver for Handpiece, medium
066.4102		Implant driver for Handpiece, long

TLX BLX A Module Set-up: Tissue Level X, Bone Level X


A Module
041.761




TLX BLX B Module Set-up: Tissue Level X, Bone Level X

Fully Tapered Tray 041.777


Implant Driver for Ratchet
066.4201
066.4207
066.4202




Implant Driver for Handpiece
066.4101
066.4107
066.4102




Alignment Pin
046.799




Depth Gauge
046.801




Depth Gauge
046.803




Depth Gauge
046.805




Depth Gauge
046.807




Depth Gauge
046.800




Depth Gauge
046.802




Depth Gauge
046.804




Depth Gauge
046.806






	Ø 3.5	Ø 3.75	Ø 4.0	Ø 4.5	Ø 5.0	Ø 5.5	Ø 6.5
soft	● 2.2	● 2.8	● 2.8	● 2.8	○ 3.2	○ 3.2	● 3.5
medium	● 2.8 (● 3.5)	○ 3.2 (● 3.7)	● 3.5 (● 3.7)	● 3.7	● 4.2 (● 4.7)	● 4.7	● 5.2 (● 6.2)
hard	○ 3.2 (● 3.5)	● 3.5 (● 3.7)	● 3.7	● 4.2	● 4.7	● 5.2	● 6.2


Round Bur
044.004
044.003




Needle Drill
026.0056




X VeloDrill™
066.1701
066.1301




X VeloDrill™
066.1702
066.1302




X VeloDrill™
066.1703
066.1303




X VeloDrill™
066.1704
066.1304




X VeloDrill™
066.1705
066.1305




X VeloDrill™
066.1706
066.1306




X VeloDrill™
066.1707
066.1307



X VeloDrill™
066.1308



X VeloDrill™
066.1309





Order List: Tissue Level X Guided, Bone Level X Guided

A MODULE – ORDER LIST

TLX/BLX Guided		
Art. No.	Image	Product
041.761		Straumann® Modular Cassette, A Module
041.766		A Modul Ratchet Tray
046.119		Ratchet
066.1100		Torque Control Device
046.064		Holding Key
041.764		Grommet Tray, 3 small + 3 large
046.400		SCS Screwdriver for Ratchet, extra-short
046.401		SCS Screwdriver for Ratchet, short
046.402		SCS Screwdriver for Ratchet, long
046.410		SCS Screwdriver for Handpiece, extra-short
046.411		SCS Screwdriver for Handpiece, short
046.412		SCS Screwdriver for Handpiece, long
041.771		Implant Depth Gauge Tray
066.2000		Implant Depth Gauge
041.762		Grommet Tray, 6 small
026.0022		BL Bone Profiler 1
026.0023		BL Bone Profiler 2
026.0024		BL Bone Profiler 3
044.304		Bone Profiler RN/RT

044.305		Bone Profiler WN/WT
041.762		A Module, Grommet Tray, 6 small
036.3300		RT Profile Drill, short, for implants Ø3.75/Ø4.8 mm
036.3301		RT Profile Drill, long, for implants Ø3.75/Ø4.8 mm
036.3302		WT Profile Drill, short, for implants Ø5.5/Ø6.5 mm
036.3303		WT Profile Drill, long, for implants Ø5.5/Ø6.5 mm

For additional tools, tweezers, AS Screwdriver refer to page 3 *A Module full list*.

B MODULE – ORDER LIST

TLX/BLX Guided		
Art. No.	Image	Product
041.776		Straumann® Modular Cassette, B Module, Base + Lid
041.780		B Module, Fully Tapered Guided Tray
034.010		Mucosa Punch, Ø3.4 mm
034.011		Mucosa Punch, Ø4.0 mm
034.012		Mucosa Punch, Ø4.7 mm
034.215		Milling Cutter, Ø2.8 mm
034.415		Milling Cutter, Ø3.5 mm
034.615		Milling Cutter, Ø4.2 mm
066.1301		X VeloDrill™, short, Ø2.2 mm
066.1302		X VeloDrill™, short, Ø2.8 mm
066.1303		X VeloDrill™, short, Ø3.2 mm
066.1304		X VeloDrill™, short, Ø3.5 mm
066.1305		X VeloDrill™, short, Ø3.7 mm
066.1306		X VeloDrill™, short, Ø4.2 mm
066.1307		X VeloDrill™, short, Ø4.7 mm



Order List: Tissue Level X Guided, Bone Level X Guided

066.1308		X VeloDrill™, short, Ø5.2 mm
066.1309		X VeloDrill™, short, Ø6.2 mm
066.1501		X VeloDrill™, medium, Ø2.2 mm
066.1502		X VeloDrill™, medium, Ø2.8 mm
066.1503		X VeloDrill™, medium, Ø3.2 mm
066.1504		X VeloDrill™, medium, Ø3.5 mm
066.1505		X VeloDrill™, medium, Ø3.7 mm
066.1506		X VeloDrill™, medium, Ø4.2 mm
066.1701		X VeloDrill™, long, Ø2.2 mm
066.1702		X VeloDrill™, long, Ø2.8 mm
066.1703		X VeloDrill™, long, Ø3.2 mm
066.1704		X VeloDrill™, long, Ø3.5 mm
066.1705		X VeloDrill™, long, Ø3.7 mm
066.1706		X VeloDrill™, long, Ø4.2 mm
066.1707		X VeloDrill™, long, Ø4.7 mm
046.799		Alignment Pin, Ø2.2 mm
046.800		Depth Gauge, Ø2.8 mm
046.801		Depth Gauge, Ø3.2 mm
046.802		Depth Gauge, Ø3.5 mm
046.803		Depth Gauge, Ø3.7 mm
046.804		Depth Gauge, Ø4.2 mm
046.805		Depth Gauge, Ø4.7 mm
046.806		Depth Gauge, Ø5.2 mm
046.807		Depth Gauge, Ø6.2 mm
066.4403		BLX Guided implant Driver for Handpiece
066.4404		BLX Guided implant Driver for Ratchet
037.3000		TLX Guided Implant Driver for ratchet, S
037.3001		TLX Guided Implant Driver for handpiece, S
037.3002		TLX Guided Implant Driver for ratchet, SP
037.3003		TLX Guided Implant Driver for handpiece, SP

C MODULE – ORDER LIST

TLX/BLX Guided		
Art. No.	Image	Product
041.772		Straumann® Modular Cassette, C Module, Guided surgery
026.0147		Drill handle, Ø2.2 mm, 1 mm/3 mm*
026.0148		Drill handle, Ø2.8 mm, 1 mm/3 mm*
026.0149		Drill handle, Ø3.2 mm, 1 mm/3 mm*
026.0150		Drill handle, Ø3.5 mm, 1 mm/3 mm*
026.0151		Drill handle, Ø3.7 mm, 1 mm/3 mm*
026.0152		Drill handle, Ø4.2 mm, 1 mm/3 mm*
034.284		Drill for Template Fixation Pin, Ø1.3 mm
034.282		Template Fixation Pin, Ø1.3 mm
034.298		Template Fixation Pin, Ø2.8/2.8 mm
034.285		Template Fixation Pin, Ø5/2.8 mm
034.286		Template Fixation Pin, Ø5/3.2 mm
034.287		Template Fixation Pin, Ø5/3.5 mm
034.288		Template Fixation Pin, Ø5/3.7 mm
034.289		Template Fixation Pin, Ø5/4.2 mm

*Drill Handles compatible with metal sleeves, for the use with PEEK sleeves please use the self-locking handles with article numbers 034.291 - 034.296.



A Module Set-up: Tissue Level X Guided, Bone Level X Guided

A Module
041.761

Ratchet Tray
041.766


Ratchet
046.119

Torque
Control
066.1100

Holding Key
046.064

Implant Depth Gauge Tray
041.771

Implant Depth Gauge
066.2000



SCS Screwdriver
for ratchet
046.400
046.401
046.402

SCS Screwdriver
for Handpiece
046.410
046.411
046.412

SCS Screwdriver Set
Grommet Tray 3 small + 3 large
041.764

Bone Profiler
026.0022
026.0023
026.0024
044.304
044.305

Bone Profiler Set
Grommet Tray 6 small
041.762




Profile drill
036.3300
036.3301
036.3302
036.3303


RT/WT Profile drill set
Grommet Tray, 6 small
041.762



B Module Set-up: Tissue Level X Guided, Bone Level X Guided


Fully Tapered Guided Tray 041.780





C Module Set-up: Tissue Level X Guided, Bone Level X Guided

C Module Guided Surgery
041.772



NOTES

Lined area for notes



International Headquarters

Institut Straumann AG
Peter Merian-Weg 12
CH-4002 Basel, Switzerland
Phone +41 (0)61 965 11 11
Fax +41 (0)61 965 11 01
www.straumann.com

© Institut Straumann AG, 2023. All rights reserved.
Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

